

**GENERAL NOTES:**

**NOTE:** The purpose of this layout drawing is primarily/specifically to illustrate the location(s) of the Siemens-supplied equipment. All non-Siemens and/or future and/or existing equipment (e.g. millwork, furniture, carts, etc.) is shown for representational purposes and may not accurately represent the actual on-site configurations. It is the responsibility of the customer/contractor that all related codes, policies (e.g. hospital requirements, building codes) and clearances to the Siemens equipment are followed when locating these ancillary items.

**NOTE:** All construction room upgrades (such as: new walls, doors, windows, millwork, plumbing, furniture, medical equip, and etc.) are shown for proposal purposes only. The customer is to provide these room upgrades. In the case of a turnkey, these room upgrades will be detailed in the contractor's quote. The contractor's quote takes precedence over the room layout drawings.

**NOTE:** This layout is preliminary only. Siemens reserves the right to make changes as dictated by new technical developments.

**NOTE:** The customer is responsible to have the transport route for Siemens equipment verified by the Structural Engineer (from point of delivery to final location).

**NOTE:** Additional anchoring may be required if the site is deemed as 'Seismic Zone' by local codes. The customer must verify the local code requirements and advise the Siemens Project Manager.

**NOTE:** Siemens is not an Architectural or Engineering firm. Drawings supplied by Siemens are not construction drawings. Therefore, these drawings are to be used only for information to complement actual construction drawings available from a customer-appointed Architectural representative or a Customer's Engineering design group. The Customer's Architect and General Contractor shall be ultimately responsible for compliance with all applicable codes and professional design requirements.

**NOTE:** Additional accessories that typically accompany the system will require storage space.

**NOTE:** All drains will need to be confirmed as meeting the requirement for drain flow and capacity as determined in the Siemens specifications by the customer.


**NOTE:** Customer water supply must meet Siemens specifications prior to instrument installation for both water specific requirements as well as meeting instrument demand specifications.

**Responsibilities:**

The Customer is responsible for site preparation for the installation of the Automation system and accessory UPS, Transformer, and Compressor, provided by Siemens. The customer also installs the UPS Transformer, and Compressor, along with all wiring breakers panels and switches for these accessory items following local code and state regulations. The Customer is responsible for the installation of the quantity and required power water and air in designated locations. Power poles are typical but not a specific requirement. The Customer also provides for the installation of Cat 5 Data lines in the locations called out in the IT diagram.

Electricians are responsible for following all State and local building codes and National Electrical codes.

**NOTE:** Contractor or Siemens to site verify locations of all existing walls. The locations of existing walls may be different than as shown in this layout proposal.



**Utility Locations 1 & 2**

Atellica System Power (2 Locations)  
 3- Total 208VAC per location  
 1- Total 120VAC per location

Atellica System Power With Eaton 9PX-3000 GL RT UPS

Atellica Sampler Handler Power  
 1- 200-240VAC 60hz 1-Phase  
 20 Amp Dedicated / Non Isolated Circuit  
 Receptacle NEMA L6-20R

Accessory Power  
 1- 120 VAC 60 Hz 1-phase  
 20 Amp Circuit  
 Receptacle Nema #5-20R

Atellica Chemistry Analyzer Power  
 1- 200-240VAC 60hz 1-Phase  
 20 Amp Dedicated / Non Isolated Circuit  
 Receptacle NEMA L6-20R

Atellica Immunoassy Analyzer Power  
 1- 200-240VAC 60hz 1-Phase  
 20 Amp Dedicated / Non Isolated Circuit  
 Receptacle NEMA L6-20R

DI Water in this location  
 Reagent Grade DI Water (CLSI C3-A4) required. Regulated to 5-30PSI with two shutoff valves, and two 3/8" compression fittings to supply the analytical modules.

Data  
 1 Laboratory Network Connection Required

**Utility Location 3**

Work Station Power:  
 3 Circuits 120VAC 20Amp  
 Duplex Receptacles  
 NEMA #5-20R

Data  
 3 - Total 3 Laboratory Network Connections

6- Total Cat 5 / 6 Point-to-Point connections to utility location 6.

**Utility Location 4**

Aptio System Power  
 From UPS, Transformer,  
 and Cut-off Switch

**Utility Location 5**

Compressed Air  
 From provided Air Compressor  
 3/4" Tubing Run From Source  
 3/4" Ball Valve with 1/2" female threaded fitting and Pressure Gauge 150 PSI min  
 House Air Must Be Clean and Oil Free

**Utility Location 6**

Accessory Power  
 120VAC 20 Amp Circuit  
 Receptacle Quad Nema 5-20R

Data Lines from Network and Workstations

No	REVISIONS	CHK'D

**SUBMITTED FOR:**

Approval  
 Your use  
 Your records  
 Installation

APPROVED: \_\_\_\_\_

**NOTES:**  
 ALL DIMENSIONS SHOWN MUST BE VERIFIED ON THE SITE. THIS DRAWING IS NOT TO BE REPRODUCED WITHOUT WRITTEN AUTHORIZATION.

**CUSTOMER:**  
 St. Luke's Hospital- East  
 100 NE St. Lukes Blvd.  
 Lee's Summit, MO  
 64086

**EQUIPMENT:** Aptio Automation  
**PG#:** 1 of 5  
**DRAWN BY:** CMH **SCALE:** 1:35  
**CHK'D BY:** \_\_\_\_\_ **DATE:** 09/19/2019  
**HOSPITAL #:** \_\_\_\_\_ **PROJECT #:** \_\_\_\_\_  
**DRAWING/SKETCH #:** \_\_\_\_\_ **REVISION #:** \_\_\_\_\_

IN ORDER TO AVOID DELAY IN INSTALLATION, SIEMENS PROJECT MANAGER SHOULD BE CONSULTED PRIOR TO INSTALLATION. FINAL ARCHITECTURAL DWGS SHOULD BE MADE AVAILABLE TO SIEMENS PROJECT MANAGER AT THIS TIME TO VERIFY THAT ALL REQUIREMENTS HAVE BEEN ADHERED TO.

THIS DRAWING DOES NOT PROVIDE DETAILS ON THE ORIGIN OF THE ELECTRICAL FEEDS FOR SIEMENS EQUIPMENT, ONLY THE CONNECTION/TERMINATION POINT. THE ORIGIN OF THE FEED WILL NEED TO BE DETERMINED BY THE CUSTOMER AND MAY REQUIRE ADDITIONAL ELECTRICAL CONSULTATION.

ANY CURRENT EQUIPMENT DISPLACED BY THE INSTALLATION OF NEW SIEMENS EQUIPMENT WILL BE THE RESPONSIBILITY OF THE CUSTOMER TO DETERMINE NEW LOCATION AND THE COORDINATION OR MOVE OF ANY REQUIRED SERVICES TO KEEP EQUIPMENT OPERATIONAL. LOCATIONS CAN BE DETERMINED IN CONSULTATION WITH SIEMENS PROJECT MANAGER.

METRIC to IMPERIAL CONVERSIONS:  
 1000mm = 39.37" 1'-0" = 304.8mm 1kg = 2.205lbs.  
 Scale 1:100

ALL DIMENSIONS SHOWN ARE IN MILLIMETERS, UNLESS OTHERWISE SPECIFIED.  
 ALL DIMENSIONS SHOWN ARE FROM FINISHED SURFACES.

THIS PLANNING PROPOSAL, TOGETHER WITH ANY ENCLOSED DOCUMENTATION AND SPECIFICATIONS, IS BASED ON THE MOST CURRENT TECHNICAL INFORMATION AVAILABLE AT THE TIME OF ISSUE. WE RESERVE THE RIGHT TO MAKE CHANGES AS DICTATED BY TECHNICAL DEVELOPMENTS.

SEISMIC REQUIREMENTS:  
 ADDITIONAL ANCHORING MAY BE REQUIRED IF SITE IS DEEMED AS 'SEISMIC ZONE' BY LOCAL CODES. THE CUSTOMER MUST VERIFY THE LOCAL CODES REQUIREMENTS AND ADVISE THE SIEMENS PROJECT MANAGER.

**GENERAL NOTES:**

**NOTE:** The purpose of this layout drawing is primarily/specifically to illustrate the location(s) of the Siemens-supplied equipment. All non-Siemens and/or future and/or existing equipment (e.g. millwork, furniture, carts, etc.) is shown for representational purposes and may not accurately represent the actual on-site configurations. It is the responsibility of the customer/contractor that all related codes, policies (e.g. hospital requirements, building codes) and clearances to the Siemens equipment are followed when locating these ancillary items.

**NOTE:** All construction room upgrades (such as: new walls, doors, windows, millwork, plumbing, furniture, medical equip, and etc.) are shown for proposal purposes only. The customer is to provide these room upgrades. In the case of a turnkey, these room upgrades will be detailed in the contractor's quote. The contractor's quote takes precedence over the room layout drawings.

**NOTE:** This layout is preliminary only. Siemens reserves the right to make changes as dictated by new technical developments.

**NOTE:** The customer is responsible to have the transport route for Siemens equipment verified by the Structural Engineer (from point of delivery to final location).

**NOTE:** Additional anchoring may be required if the site is deemed as 'Seismic Zone' by local codes, the customer must verify the local code requirements and advise the Siemens Project Manager.

**NOTE:** Siemens is not an Architectural or Engineering firm. Drawings supplied by Siemens are not construction drawings. Therefore, these drawings are to be used only for information to complement actual construction drawings available from a customer-appointed Architectural representative or a Customer's Engineering design group. The Customer's Architect and General Contractor shall be ultimately responsible for compliance with all applicable codes and professional design requirements.

**NOTE:** Additional accessories that typically accompany the system will require storage space.

**NOTE:** All drains will need to be confirmed as meeting the requirement for drain flow and capacity as determined in the Siemens specifications by the customer.

**NOTE:** Customer water supply must meet Siemens specifications prior to instrument installation for both water specific requirements as well as meeting instrument demand specifications.

**Responsibilities:**

The Customer is responsible for site preparation for the installation of the Automation system and accessory UPS, Transformer, and Compressor, provided by Siemens. The customer also installs the UPS Transformer, and Compressor, along with all wiring breakers panels and switches for these accessory items following local code and state regulations.


The Customer is responsible for the installation of the quantity and required power water and air in designated locations.

Power poles are typical but not a specific requirement.

The Customer also provides for the installation of Cat 5 Data lines in the locations called out in the IT diagram.

Electricians are responsible for following all State and local building codes and National Electrical codes.

**NOTE:** Contractor or Siemens to site verify locations of all existing walls. The locations of existing walls may be different than as shown in this layout proposal.



No	REVISIONS	CHKD

**SUBMITTED FOR:**

Approval

Your use

Your records

Installation

**APPROVED:** \_\_\_\_\_

**NOTES:**  
ALL DIMENSIONS SHOWN MUST BE VERIFIED ON THE SITE. THIS DRAWING IS NOT TO BE REPRODUCED WITHOUT WRITTEN AUTHORIZATION.

**CUSTOMER:**  
St. Lukes Hospital- East  
100 NE St. Lukes Blvd  
Lee's Summit, MO  
64086

<b>EQUIPMENT:</b>	Aptio Automation
<b>PG#:</b>	2 of 5
<b>DRAWN BY:</b>	CMH
<b>CHKD BY:</b>	DATE: 09/19/2019
<b>HOSPITAL #</b>	<b>PROJECT #</b>

IN ORDER TO AVOID DELAY IN INSTALLATION, SIEMENS PROJECT MANAGER SHOULD BE CONSULTED PRIOR TO INSTALLATION. FINAL ARCHITECTURAL DWGS SHOULD BE MADE AVAILABLE TO SIEMENS PROJECT MANAGER AT THIS TIME TO VERIFY THAT ALL REQUIREMENTS HAVE BEEN ADHERED TO.

THIS DRAWING DOES NOT PROVIDE DETAILS ON THE ORIGIN OF THE ELECTRICAL FEEDS FOR SIEMENS EQUIPMENT, ONLY THE CONNECTION/TERMINATION POINT. THE ORIGIN OF THE FEED WILL NEED TO BE DETERMINED BY THE CUSTOMER AND MAY REQUIRE ADDITIONAL ELECTRICAL CONSULTATION.

ANY CURRENT EQUIPMENT DISPLACED BY THE INSTALLATION OF NEW SIEMENS EQUIPMENT WILL BE THE RESPONSIBILITY OF THE CUSTOMER TO DETERMINE NEW LOCATION AND THE COORDINATION OR MOVE OF ANY REQUIRED SERVICES TO KEEP EQUIPMENT OPERATIONAL. LOCATIONS CAN BE DETERMINED IN CONSULTATION WITH SIEMENS PROJECT MANAGER.

**METRIC to IMPERIAL CONVERSIONS:**  
1000mm = 39.37" 1'-0" = 304.8mm 1kg = 2.205lbs.  
Scale 1:100

ALL DIMENSIONS SHOWN ARE IN MILLIMETERS, UNLESS OTHERWISE SPECIFIED.  
ALL DIMENSIONS SHOWN ARE FROM FINISHED SURFACES.

THIS PLANNING PROPOSAL, TOGETHER WITH ANY ENCLOSED DOCUMENTATION AND SPECIFICATIONS, IS BASED ON THE MOST CURRENT TECHNICAL INFORMATION AVAILABLE AT THE TIME OF ISSUE. WE RESERVE THE RIGHT TO MAKE CHANGES AS DICTATED BY TECHNICAL DEVELOPMENTS.

**SEISMIC REQUIREMENTS:**  
ADDITIONAL ANCHORING MAY BE REQUIRED IF SITE IS DEEMED AS 'SEISMIC ZONE' BY LOCAL CODES. THE CUSTOMER MUST VERIFY THE LOCAL CODES REQUIREMENTS AND ADVISE THE SIEMENS PROJECT MANAGER.

**NOTE:** The purpose of this layout drawing is primarily to illustrate the location(s) of the Siemens-supplied equipment. All non-Siemens and/or future and/or existing equipment (e.g. millwork, furniture, carts, etc.) is shown for representational purposes and may not accurately represent the actual on-site configurations. It is the responsibility of the customer/contractor that all related codes, policies (e.g. hospital requirements, building codes) and clearances to the Siemens equipment are followed when locating these ancillary items.

**NOTE:** All construction room upgrades (such as: new walls, doors, windows, millwork, plumbing, furniture, medical equip, and etc.) are shown for proposal purposes only. The customer is to provide these room upgrades. In the case of a turnkey, these room upgrades will be detailed in the contractor's quote. The contractor's quote takes precedence over the room layout drawings.

**NOTE:** This layout is preliminary only. Siemens reserves the right to make changes as dictated by new technical developments.

**NOTE:** The customer is responsible to have the transport route for Siemens equipment verified by the Structural Engineer (from point of delivery to final location).

**NOTE:** Additional anchoring may be required if the site is deemed as "Seismic Zone" by local codes. The customer must verify the local code requirements and advise the Siemens Project Manager.

**NOTE:** Siemens is not an Architectural or Engineering firm. Drawings supplied by Siemens are not construction drawings. Therefore, these drawings are to be used only for information to complement actual construction drawings available from a customer-appointed Architectural representative or a Customer's Engineering design group. The Customer's Architect and General Contractor shall be ultimately responsible for compliance with all applicable codes and professional design requirements.

**NOTE:** Additional accessories that typically accompany the system will require storage space.

**NOTE:** All drains will need to be confirmed as meeting the requirement for drain flow and capacity as determined in the Siemens specifications by the customer.

**NOTE:** Customer water supply must meet Siemens specifications prior to instrument installation for both water specific requirements as well as meeting instrument demand specifications.

**Responsibilities:**  
The Customer is responsible for site preparation for the installation of the Automation system and accessory UPS, Transformer, and Compressor, provided by Siemens. The customer also installs the UPS Transformer, and Compressor, along with all wiring breakers panels and switches for these accessory items following local code and state regulations. The Customer is responsible for the installation of the quantity and required power water and air in designated locations. Power poles are typical but not a specific requirement. The Customer also provides for the installation of Cat 5 Data lines in the locations called out in the IT diagram. Electricians are responsible for following all State and local building codes and National Electrical codes.

**Unrestricted**

**NOTE:** Contractor or Siemens to site verify locations of all existing walls. The locations of existing walls may be different than as shown in this layout proposal.

System Utilities					
Category	Title	Requirement	Qty	Utility Location	Comments
Power	SF4 Compressor Main Power	208/230 VAC 60 Hz, 3-phase dedicated circuit, 30 Amp circuit with 25/25 amp fast acting fuse in the cutoff switch. Respectably	1	TBD	AWG 10 Wire Dedicated Circuit Emergency Power
Power	Aptio UPS Main Power (1-phase 12kVa)	4-wire 208 VAC 60Hz with neutral, Single Phase, 100 AMP 2 AWG	1	TBD	Dedicated Circuit Emergency Power Neutral wire required.
Power	Power connections from UPS to Track	1-P (3-wire, L1, L2, G) See block diagram page	1	4	Connections include transformer, Shut Off Switch inside Lab, and Aptio Track
Power	Command Center Power	120VAC 1-phase 20A, Duplex NEMA #5-20R	3	3	Dedicated Circuit Emergency Power
Power	Aptio Auxiliary power	120VAC 1-phase 20A, Duplex NEMA #5-20R	1	6	Dedicated Circuit Emergency Power
Power	Atellica Sample Handler	208VAC 1-phase 20A Long Delay Breaker, HBL5461 or or equivalent NEMA 6-20R	2	1 & 2	Refer to Atellica Site Survey for more detail
Power	Atellica Auxiliary power	120VAC 1-phase 20A, Duplex NEMA #5-20R	2	1 & 2	Refer to Atellica Site Survey for more detail
Power	Atellica Chemistry Analyzer	208VAC 1-phase 20A Long Delay Breaker, HBL5461 or or equivalent NEMA 6-20R	2	1 & 2	Refer to Atellica Site Survey for more detail
Power	Atellica Immunoassay Analyzer	208VAC 1-phase 20A Long Delay Breaker, HBL5461 or or equivalent NEMA 6-20R	2	1 & 2	Refer to Atellica Site Survey for more detail
Data	Command Center Data	RJ-45 Connections to Laboratory Network	3	3	Includes connection for Hospital LIS Computer
Data	Command Center Point-to-Point	RJ-45 Connections from Command Center to Aptio Track	6	3 & 6	Connections from Command Center to Aptio Track
Water	Atellica Chemistry Analyzer	Reagent Grade Water (CLSI C3-A4) required in this location regulated to 5-30 psi with shutoff valves, and 3/8" compression fittings to supply the analytical modules.	2	1 & 2	Refer to Atellica Site Survey for more detail
Water	Atellica Immunoassay Analyzer	Reagent Grade Water (CLSI C3-A4) required in this location regulated to 5-30 psi with shutoff valves, and 3/8" compression fittings to supply the analytical modules.	2	1 & 2	Refer to Atellica Site Survey for more detail
Air	Copper tubing routed from the Compressor to the Aptio track	Air Supply Line copper tubing run with 3/4 inch female threaded fitting at end near the track.	1	5	3/4 inch copper tubing routed from the Compressor to the Aptio track.
Air	Supply Shutoff valve, air gauge 150 PSI min, T-fitting, NPT 3/4 inch female fitting, and drain valve.		1	5	Located in the LAB next to the Aptio Automation. Refer to pole location.

### Aptio Automation Environmental Impact

INPECO Technical Specifications Vers. 1.38 LDAT-000.835.02.01.02(002)

		Metric Units		Standard Units
Weight of Aptio modules	2640	Kg	5820	Lbs
Weight of separate analyzer	3296	Kg	7266	Lbs
Aptio Compressed air (NL/min)	211	NL/min		
Total power consumption of modules through internal busway (VA)	8	KVA		
Total power consumption (VA) of separate analyzers**	12	VA		
Heat of modules with analyzers (BTU/h)	53738	BTU/h		
Waste drain capacity	110.0	L/h	29.1	G/h
DI water require-ments	78.0	L/h	20.6	G/h

The calculations here are based on the site surveys for the equipment shown at publication and are subject to change.

\*\* Analyzers are not included in the total busway calculation. Power is nominal and does not include Inrush Current.

#### Floor Loading Capacity Requirements:


Unless otherwise noted, all floor loading capacity requirements are for less than 300 Kg/square meter.

#### Environmental:

Temperature: from 5°C to 40°C.  
Humidity: Max: 80% at 31°C to 50% at 40°C.

#### Compressed Air:

Supply air pressure range: 7 bar – 9 bar.  
Air quality: Output air quality for all family members will meet ISO8573-1:2010 [1:4:1] (clean, oil-free, and dry compressed air).  
ISO 8573 states:  
Less than 20,000 particles per cubic meter for particle sizes between 0.1 micrometer and 0.5 micrometer  
Less than 400 particles per cubic meter for particle sizes between 0.5 micrometer and 1.0 micrometer  
Less than 10 particles per cubic meter for particle sizes between 1.0 micrometer and 5.0 micrometer  
Pressure dewpoint less than or equal to +3 °C  
Total oil less than or equal to 0.01 milligrams per cubic meter  
Total compressed air is adjusted by a factor of 1.5.  
Reference conditions for Normal Liters per minute (NL/min) in this spreadsheet: 20 degrees C temperature and 1 bar pressure (per ISO 1217).



A = 6' B = 5' 10" C = 3'



SUBMITTED FOR:

Approval  
 Your use  
 Your records  
 Installation

APPROVED: \_\_\_\_\_

NOTES:  
ALL DIMENSIONS SHOWN MUST BE VERIFIED ON THE SITE. THIS DRAWING IS NOT TO BE REPRODUCED WITHOUT WRITTEN AUTHORIZATION.

CUSTOMER:  
**St. Lukes Hospital- East  
100 NE St. Lukes Blvd  
Lee's Summit, MO  
64086**

EQUIPMENT: **Aptio Automation**  
PAGE: **3 of 5**  
DRAWN BY: **CMH** SCALE:  
CHKD BY: DATE: **09/19/2019**

**GENERAL NOTES:**

**NOTE:** The purpose of this layout drawing is primarily to illustrate the location(s) of the Siemens-supplied equipment. All non-Siemens and/or future and/or existing equipment (e.g. millwork, furniture, carts, etc.) is shown for representational purposes and may not accurately represent the actual on-site configurations. It is the responsibility of the customer/contractor that all related codes, policies (e.g. hospital requirements, building codes) and clearances to the Siemens equipment are followed when locating these ancillary items.

**NOTE:** All construction room upgrades (such as: new walls, doors, windows, millwork, plumbing, furniture, medical equip, and etc.) are shown for proposal purposes only. The customer is to provide these room upgrades. In the case of a turnkey, these room upgrades will be detailed in the contractor's quote. The contractor's quote takes precedence over the room layout drawings.

**NOTE:** This layout is preliminary only. Siemens reserves the right to make changes as dictated by new technical developments.

**NOTE:** The customer is responsible to have the transport route for Siemens equipment verified by the Structural Engineer (from point of delivery to final location).

**NOTE:** Additional anchoring may be required if the site is deemed as 'Seismic Zone' by local codes, the customer must verify the local code requirements and advise the Siemens Project Manager.

**NOTE:** Siemens is not an Architectural or Engineering firm. Drawings supplied by Siemens are not construction drawings. Therefore, these drawings are to be used only for information to complement actual construction drawings available from a customer-appointed Architectural representative or a Customer's Engineering design group. The Customer's Architect and General Contractor shall be ultimately responsible for compliance with all applicable codes and professional design requirements.

**NOTE:** Additional accessories that typically accompany the system will require storage space.

**NOTE:** All drains will need to be confirmed as meeting the requirement for drain flow and capacity as determined in the Siemens specifications by the customer.

**NOTE:** Customer water supply must meet Siemens specifications prior to instrument installation for both water specific requirements as well as meeting instrument demand specifications.

**Responsibilities:**  
The Customer is responsible for site preparation for the installation of the Automation system and accessory UPS, Transformer, and Compressor, provided by Siemens.

The customer also installs the UPS Transformer, and Compressor, along with all wiring breakers panels and switches for these accessory items following local code and state regulations.

The Customer is responsible for the installation of the quantity and required power water and air in designated locations. Power poles are typical but not a specific requirement.

The Customer also provides for the installation of Cat 5 Data lines in the locations called out in the IT diagram.

Electricians are responsible for following all State and local building codes and National Electrical codes

**Aptio Main Power**




**40 Amp breaker or 40 Amp Disconnect Switch**



**Eaton 9155 UPS**



**Step-Up Transformer**



Internal circuit breakers should not be used for customer main disconnect.


**NOTES**  
All connections are hardwired.  
Installation must conform to local electrical codes.  
Loads should be considered non-linear.

UPS wiring must be copper rated at 600 VAC, 75° C  
UPS wiring recommendations are extracted from the Eaton User Manuals referenced below.

All voltages are listed as nominal voltages.  
Allowable fluctuations = +/- 10%

Refer to the Aptio Environmental Spreadsheet (DCAU-F00.811.02, Compressor and UPS Sizing tab) for UPS sizing and associated SMNs.


Refer also to the latest version of the Eaton User Manual:  
<http://powerquality.eaton.com/Products-services/Backup-Power-UPS/9155.aspx#documentationtab>



**Figure 18. UPS with MBM Rear View**

Aptio™ Automation is a trademark of Siemens Healthcare Diagnostics.


Eaton 9155 1-phase 12 kVA UPS and 1-Phase Transformer		
9155 UPS Weight and Dimensions	9155 UPS Electrical Requirements	1-Phase Transformer
<b>Weight:</b> 171 (kg) 377 (lbs)	4-wire 1-phase dedicated circuit direct wired to: • 12 kVA UPS (neutral required at the UPS) • step-up transformer • Electrical disconnect switch or 40 Amp breaker • Aptio power plate	<b>Weight:</b> 17 (kg) 38 (lbs)
<b>Dimensions:</b> (H x W x D) 32" x 12" x 33"	<b>Voltage:</b> 208VAC 60Hz	<b>Dimensions:</b> (H x W x D) 13.5" x 5.5" x 5.13"
	<b>Current:</b> 100 Amp Circuit	<b>Heat Output:</b> 812 BTU's at full capacity
<b>Heat Output:</b> 4,512 BTU/Hr	Hard Wired 2 AWG.	



**Power Protection & Power Quality Specialists**

Primary configuration: 1 winding – 4 sections  
Single winding where the voltage delta is dependent on where the HV & LV conductors are connected

Wiring Diagram – 208VAC Input (Wiring Instructions to be provided with transformer when shipped):



- Wiring Instructions:**
1. Connect H4 together with one of your input wires (lower voltage) and one of your output wires (higher voltage). You should have a total of 3 wires combined together for this one connection.
  2. Connect H2 and H3 together. You should have a total of 2 wires combined together for this one connection.
  3. Connect your second input wire (lower voltage) to H1 and X4. You should have a total of 3 wires combined together for this one connection.
  4. Connect X2 and X3 together. You should have a total of 2 wires combined together for this one connection.
  5. Connect your second output wire (higher voltage) together with X1. You should have a total of 2 wires combined together for this one connection.
  6. Re-verify wires in previous steps. Check that all connections are tight and well insulated.
  7. Prior to powering the Siemens Aptio, apply power to the transformer and verify that you have the proper voltage (229VAC) for the Siemens Aptio.

Secondary configuration: 220VAC Input – see below

Winding material: Copper

Temperature rise: 115° C

Electrostatically shielded: No

1590 North Roberts Road - Suite 104 Kennesaw, Georgia 30144 (770) 499-0900 – Phone (770) 499-9493 – Fax [www.power-place.com](http://www.power-place.com) – Website

Page 2      Rev 1 – 01/16/14  
Distributors of: - UPS Systems  
- Power Line Monitor Services  
- Surge Protectors  
- Power Conditioners  
- Generator Systems

No REVISIONS CHKD

SUBMITTED FOR:  
 Approval  
 Your use  
 Your records  
 Installation  
APPROVED: \_\_\_\_\_



**NOTES:**  
ALL DIMENSIONS SHOWN MUST BE VERIFIED ON THE SITE. THIS DRAWING IS NOT TO BE REPRODUCED WITHOUT WRITTEN AUTHORIZATION.

EQUIPMENT: Eaton 9155 12 kVa UPS	
FOR: Page 4	
DRAWN BY: KM	SCALE:
CHKD BY: BT	DATE: 05/02/2019
HOSPITAL #	PROJECT #

IN ORDER TO AVOID DELAY IN INSTALLATION, SIEMENS PROJECT MANAGER SHOULD BE CONSULTED PRIOR TO INSTALLATION. FINAL ARCHITECTURAL DWGS SHOULD BE MADE AVAILABLE TO SIEMENS PROJECT MANAGER AT THIS TIME TO VERIFY THAT ALL REQUIREMENTS HAVE BEEN ADHERED TO.

THIS DRAWING DOES NOT PROVIDE DETAILS ON THE ORIGIN OF THE ELECTRICAL FEEDS FOR SIEMENS EQUIPMENT. ONLY THE CONNECTION/TERMINATION POINT, THE ORIGIN OF THE FEED WILL NEED TO BE DETERMINED BY THE CUSTOMER AND MAY REQUIRE ADDITIONAL ELECTRICAL CONSULTATION.

ANY CURRENT EQUIPMENT DISPLACED BY THE INSTALLATION OF NEW SIEMENS EQUIPMENT WILL BE THE RESPONSIBILITY OF THE CUSTOMER TO DETERMINE NEW LOCATION AND THE COORDINATION OR MOVE OF ANY REQUIRED SERVICES TO KEEP EQUIPMENT OPERATIONAL. LOCATIONS CAN BE DETERMINED IN CONSULTATION WITH SIEMENS PROJECT MANAGER.

**METRIC TO IMPERIAL CONVERSIONS:**  
1000mm = 39.37"    1'-0" = 304.8mm    1kg = 2.205lbs.  
Scale: 1:100

ALL DIMENSIONS SHOWN ARE IN MILLIMETERS, UNLESS OTHERWISE SPECIFIED.  
ALL DIMENSIONS SHOWN ARE FROM FINISHED SURFACES.

THIS PLANNING PROPOSAL, TOGETHER WITH ANY ENCLOSED DOCUMENTATION AND SPECIFICATIONS, IS BASED ON THE MOST CURRENT TECHNICAL INFORMATION AVAILABLE AT THE TIME OF ISSUE. WE RESERVE THE RIGHT TO MAKE CHANGES AS DICTATED BY TECHNICAL DEVELOPMENTS.

**SEISMIC REQUIREMENTS:**  
ADDITIONAL ANCHORING MAY BE REQUIRED IF SITE IS DEEMED AS 'SEISMIC ZONE' BY LOCAL CODES. THE CUSTOMER MUST VERIFY THE LOCAL CODES REQUIREMENTS AND ADVISE THE SIEMENS PROJECT MANAGER.

**GENERAL NOTES:**

**NOTE:** The purpose of this layout drawing is primarily to illustrate the location(s) of the Siemens-supplied equipment. All non-Siemens and/or future and/or existing equipment (e.g. millwork, furniture, carts, etc.) is shown for representational purposes and may not accurately represent the actual on-site configurations. It is the responsibility of the customer/contractor that all related codes, policies (e.g. hospital requirements, building codes) and clearances to the Siemens equipment are followed when locating these ancillary items.

**NOTE:** All construction room upgrades (such as: new walls, doors, windows, millwork, plumbing, furniture, medical equip, and etc.) are shown for proposal purposes only. The customer is to provide these room upgrades. In the case of a turnkey, these room upgrades will be detailed in the contractor's quote. The contractor's quote takes precedence over the room layout drawings.

**NOTE:** This layout is preliminary only. Siemens reserves the right to make changes as dictated by new technical developments.

**NOTE:** The customer is responsible to have the transport route for Siemens equipment verified by the Structural Engineer (from point of delivery to final location).

**NOTE:** Additional anchoring may be required if the site is deemed as 'Seismic Zone' by local codes, the customer must verify the local code requirements and advise the Siemens Project Manager.

**NOTE:** Siemens is not an Architectural or Engineering firm. Drawings supplied by Siemens are not construction drawings. Therefore, these drawings are to be used only for information to complement actual construction drawings available from a customer-appointed Architectural representative or a Customer's Engineering design group. The Customer's Architect and General Contractor shall be ultimately responsible for compliance with all applicable codes and professional design requirements.

**NOTE:** Additional accessories that typically accompany the system will require storage space.

**NOTE:** All drains will need to be confirmed as meeting the requirement for drain flow and capacity as determined in the Siemens specifications by the customer.

**NOTE:** Customer water supply must meet Siemens specifications prior to instrument installation for both water specific requirements as well as meeting instrument demand specifications.

**Responsibilities:**  
The Customer is responsible for site preparation for the installation of the Automation system and accessory UPS, Transformer, and Compressor, provided by Siemens.




The customer also installs the UPS Transformer, and Compressor, along with all wiring breakers panels and switches for these accessory items following local code and state regulations.

The Customer is responsible for the installation of the quantity and required power water and air in designated locations. Power poles are typical but not a specific requirement.

The Customer also provides for the installation of Cat 5 Data lines in the locations called out in the IT diagram.

**NOTE:** Contractor or Siemens to site verify locations of all existing walls. The locations of existing walls may be different than as shown in this layout proposal.

Atlas Copco SF4 Compressor			
Compressor Weight and Dimensions	Room Environment Operating	Electrical Requirements	Plumbing Requirements
<b>Weight:</b> 244 (kg) 538 (lbs)	For indoor use only. Room temperature 65-85 F. <b>Adequate Room Ventilation Required</b>	1-Dedicated three phase circuit	Tubing 3/4" copper tubing run from the compressor to the Aptio Automation location.
<b>Dimensions:</b> (H x W x D) 57 3/4", 44", 25"	<b>Noise Level:</b> 58 dB (A) at 1m	<b>Voltage:</b> 208/230VAC 60Hz +/- 10%	Required accessories & fittings for air supply line.
	<b>Max Pressure:</b> 145 psi @ 10 bar	Cutoff switch with fast acting 25/25 Amp fuses.	<b>Compressor End:</b> 3/4" Shutoff Valve with 3/4" Female Threaded Fitting.
	<b>Capacity:</b> 12.5 CFM (354 NL)	<b>Current:</b> 30 Amp Circuit.	<b>Aptio End:</b> Pressure Gauge (150psi), 3/4" Female Threaded Fitting.
<b>Nominal Motor Power:</b> 5 HP	<b>Thermal Output:</b> 12000 BTU/hr Continuous	AWG 10 Hard Wired	<b>Drain:</b> Necessary at the Compressor End for Condensation.



No REVISIONS

No	REVISIONS	CHKD

SUBMITTED FOR:  
 Approval  
 Your use  
 Your records  
 Installation  
 APPROVED: \_\_\_\_\_



**NOTES:**  
 ALL DIMENSIONS SHOWN MUST BE VERIFIED ON THE SITE. THIS DRAWING IS NOT TO BE REPRODUCED WITHOUT WRITTEN AUTHORIZATION.

CUSTOMER: \_\_\_\_\_

EQUIPMENT:  
 Atlas Copco SF 4 Compressor  
 TCR: Page 5  
 DRAWN BY: KM SCALE:  
 CHKD BY: GM DATE: 03/31/2020  
 HOSPITAL # PROJECT #

IN ORDER TO AVOID DELAY IN INSTALLATION, SIEMENS PROJECT MANAGER SHOULD BE CONSULTED PRIOR TO INSTALLATION. FINAL ARCHITECTURAL DWGS SHOULD BE MADE AVAILABLE TO SIEMENS PROJECT MANAGER AT THIS TIME TO VERIFY THAT ALL REQUIREMENTS HAVE BEEN ADHERED TO.

THIS DRAWING DOES NOT PROVIDE DETAILS ON THE ORIGIN OF THE ELECTRICAL FEEDS FOR SIEMENS EQUIPMENT. ONLY THE CONNECTION/TERMINATION POINT, THE ORIGIN OF THE FEED WILL NEED TO BE DETERMINED BY THE CUSTOMER AND MAY REQUIRE ADDITIONAL ELECTRICAL CONSULTATION.

ANY CURRENT EQUIPMENT DISPLACED BY THE INSTALLATION OF NEW SIEMENS EQUIPMENT WILL BE THE RESPONSIBILITY OF THE CUSTOMER TO DETERMINE NEW LOCATION AND THE COORDINATION OR MOVE OF ANY REQUIRED SERVICES TO KEEP EQUIPMENT OPERATIONAL. LOCATIONS CAN BE DETERMINED IN CONSULTATION WITH SIEMENS PROJECT MANAGER.

METRIC to IMPERIAL CONVERSIONS:  
 1000mm = 39.37" 1'-0" = 304.8mm 1kg = 2.205lbs.  
 Scale 1:100

ALL DIMENSIONS SHOWN ARE IN MILLIMETERS, UNLESS OTHERWISE SPECIFIED.  
 ALL DIMENSIONS SHOWN ARE FROM FINISHED SURFACES.

THIS PLANNING PROPOSAL, TOGETHER WITH ANY ENCLOSED DOCUMENTATION AND SPECIFICATIONS, IS BASED ON THE MOST CURRENT TECHNICAL INFORMATION AVAILABLE AT THE TIME OF ISSUE. WE RESERVE THE RIGHT TO MAKE CHANGES AS DICTATED BY TECHNICAL DEVELOPMENTS.

SEISMIC REQUIREMENTS:  
 ADDITIONAL ANCHORING MAY BE REQUIRED IF SITE IS DEEMED AS 'SEISMIC ZONE' BY LOCAL CODES. THE CUSTOMER MUST VERIFY THE LOCAL CODES REQUIREMENTS AND ADVISE THE SIEMENS PROJECT MANAGER.

**GENERAL NOTES:**

**NOTE:** The purpose of this layout drawing is primarily/specifically to illustrate the location(s) of the Siemens-supplied equipment. All non-Siemens and/or future and/or existing equipment (e.g., millwork, furniture, carts, etc.) is shown for representational purposes and may not accurately represent the actual on-site configurations. It is the responsibility of the customer/contractor that all related codes, policies (e.g. hospital requirements, building codes) and clearances to the Siemens equipment are followed when locating these ancillary items.

**NOTE:** All construction room upgrades (such as: new walls, doors, windows, millwork, plumbing, furniture, medical equip, and etc.) are shown for proposal purposes only. The customer is to provide these room upgrades. In the case of a turnkey, these room upgrades will be detailed in the contractor's quote. The contractor's quote takes precedence over the room layout drawings.

**NOTE:** This layout is preliminary only. Siemens reserves the right to make changes as dictated by new technical developments.

**NOTE:** The customer is responsible to have the transport route for Siemens equipment verified by the Structural Engineer (from point of delivery to final location).

**NOTE:** Additional anchoring may be required if the site is deemed as "Seismic Zone" by local codes. The customer must verify the local code requirements and advise the Siemens Project Manager.

**NOTE:** Siemens is not an Architectural or Engineering firm. Drawings supplied by Siemens are not construction drawings. Therefore, these drawings are to be used only for information to complement actual construction drawings available from a customer-appointed Architectural representative or a Customer's Engineering design group. The Customer's Architect and General Contractor shall be ultimately responsible for compliance with all applicable codes and professional design requirements.

**NOTE:** Additional accessories that typically accompany the system will require storage space.

**NOTE:** All drains will need to be confirmed as meeting the requirement for drain flow and capacity as determined in the Siemens specifications by the customer.

**NOTE:** Customer water supply must meet Siemens specifications prior to instrument installation for both water specific requirements as well as meeting instrument demand specifications.


**Responsibilities:**

The Customer is responsible for site preparation for the installation of the Automation system and accessory UPS, Transformer, and Compressor, provided by Siemens. The customer also installs the UPS Transformer, and Compressor, along with all wiring breakers panels and switches for these accessory items following local code and state regulations.


The Customer is responsible for the installation of the quantity and required power water and air in designated locations. Power poles are typical but not a specific requirement. The Customer also provides for the installation of Cat 5 Data lines in the locations called out in the IT diagram.

Electricians are responsible for following all State and local building codes and National Electrical codes

**NOTE:** Contractor or Siemens to site verify locations of all existing walls. The locations of existing walls may be different than as shown in this layout proposal.



**9PX 3000GLRT**



- ① Socket for connection to AC power source
- ② Slot for optional communication card
- ③ Relay output contact
- ④ Connector for additional battery module
- ⑤ Primary group: outlets for connection of critical equipment
- ⑥ Group 1: programmable outlets for connection of equipment
- ⑦ Group 2: programmable outlets for connection of equipment
- ⑧ Connector for automatic recognition of an additional battery module
- ⑨ RS232 communication port
- ⑩ USB communication port
- ⑪ Connector for ROO (Remote On/Off) control and RPO (Remote Power Off)

**APELLICA POWER PLUG**  
EACH MODULE HAS A NEMA 6-20P PLUG STRAIGHT BLADE 20 AMP



**EATON UPS**  
UPS WITH A NEMA 6-20R AND BATTERY BACKUP 9PX3000GLRT / 9PXEBM72RT




**CUSTOMER TO INSTALL WALL PLUG**  
20 AMP NEMA L620R TURNLOCK SINGLE RECEPTACLE




H-3.5"  
W-1' 5 1/4"  
D-1' 11 3/4"

UPS/Battery Pair - 1 set per Instrument



**Rack Views**



No REVISIONS CHKD


SUBMITTED FOR:

- Approval
- Your use
- Your records
- Installation


APPROVED: \_\_\_\_\_

NOTES:  
ALL DIMENSIONS SHOWN MUST BE VERIFIED ON THE SITE. THIS DRAWING IS NOT TO BE REPRODUCED WITHOUT WRITTEN AUTHORIZATION.


CUSTOMER:	
EQUIPMENT:	Atellica Information
PAGE:	Page 6
DRAWN BY:	CMH SCALE:
CHKD BY:	DATE: 09/19/2019
HOSPITAL #:	PROJECT #:
DRAWING/SKETCH #:	REVISION #:

	NA CUSTOM WATER SYSTEMS - MILLI-Q	
	EQUIPMENT LOCATIONS AND UTILITIES	<input type="checkbox"/> APPROVED <input type="checkbox"/> REVISE AS MARKED AND RESUBMIT
SIGN: _____ COMPANY _____ DATE _____		

REVISIONS			
REV	BY	DESCRIPTION	DATE
1			



ACCESS AREA REQUIRED



ELEVATION VIEW

INTERFACE CALL-OUT TABLE			
TAG	DESCRIPTION	LOCATION	In Place
1-1	Feed water source - provide 3/4" npt male connection with isolation ball valve. 30 LPM at 50 - 85 psi. Temperature blending valve recommended if feed water temperature is expected to drop below 10°C (50 °F, max 95°F). Blending valve and valve installation by others.	On wall, within 5 feet of system. If under a sink, provide 2" Ø hole in counter or casework for tubing for connection to system by MilliporeSigma	<input type="checkbox"/> YES <input type="checkbox"/> NO
2-1	Reject from make-up system and pre-treatment	MilliporeSigma will run tubing from system to floor drain or 2" standpipe provided by customer. Standpipe max height 24 inches	<input type="checkbox"/> YES <input type="checkbox"/> NO
3-1	120VAC for water system, provide (2) 5-15R quad receptacle. Make-up system power 1,100 VA.	Locate as shown. Must be within 36" of the make-up system	<input type="checkbox"/> YES <input type="checkbox"/> NO
4-1	Distribution recirculating loop - MilliporeSigma will provide & install MAX 60 feet of 16 mm OD (10 mm ID) flexible tubing to feed analyzer(s). Conduit may be required (by others) following MilliporeSigma conduit guidelines.	Surface mounted in accessible location, conduit required if loop runs in ceiling or under the floor. Conduit required and in place?	<input type="checkbox"/> YES <input type="checkbox"/> NO
5-1	Floor drain or stand pipe by others, capable of handling 30LPM flow.	MilliporeSigma will run tubing from system to floor drain or 2" standpipe provided by customer. Standpipe max height: 24 inches	<input type="checkbox"/> YES <input type="checkbox"/> NO
6-1	Ethernet port for each system for remote access capabilities (by others)	Approx. location shown on drawing.	<input type="checkbox"/> YES <input type="checkbox"/> NO

Weight			
	Dry Weight	Shipping Weight	Operating Weight
Milli-Q™ CLX 7120	127 kg (280 lb)	155 kg (341.7 lb)	303 kg (668 lb)
Milli-Q™ CLX 7150	139 kg (306.4 lb)	143 kg (368.2 lb)	318 kg (701.1 lb)

This drawing and the data contained thereon are confidential and the property of Millipore Corp. Any other use or reproduction of this drawing is expressly prohibited without Millipore's consent.

**MILLIPORESIGMA**

NA CUSTOM WATER SYSTEMS - MILLI-Q 400 Summit Drive - Burlington MA, 01803			
CLIENT		TITLE	
Milli-Q™ CLX 7120 / 7150		Plan & Elevation View	
DESIGNED BY	DATE	SCALE	DRW NO.
XXX	XX/XX/XX	A	
CHECKED BY	DATE	SCALE	DRW NO.
XXX	XX/XX/XX	A	
REVISED BY	DATE	SCALE	DRW NO.
LOCATION	CP	STATUS	SHEET 1 OF 1

DO NOT SCALE PRINT

Atellica Solution

# General Report All Configurations

Site Survey

**Customer:**

**Address:**

**Department:**

**Room:**

**Contact Person:**

**Telephone:**

**Cust. specific no.:**

**Cust. no.:**

**Date:**



## Document Version

Siemens Healthcare GmbH reserves the right to change its products and services at any time.

In addition, manuals are subject to change without notice. The hardcopy documents correspond to the version at the time of system delivery and/or printout. Versions to hardcopy documentation are not automatically distributed.

Please contact your local Siemens Healthineers office to order a current version or refer to our website <http://www.siemens-healthineers.com>.

## Disclaimer

Siemens Healthcare GmbH provides this documentation "as is" without the assumption of any liability under any theory of law.

The content described herein shall be used by qualified personnel who are employed by Siemens Healthcare GmbH or one of its affiliates or who are otherwise authorized by Siemens Healthcare GmbH or its affiliates to use such documents.

Assemblers and other persons who are not employed by or otherwise directly affiliated with or authorized by Siemens Healthcare GmbH or one of its affiliates are not entitled to use this documentation without prior written authority.

## Copyright

"© Siemens Healthineers, 2020" refers to the copyright of a Siemens Healthineers entity such as:

Siemens Healthcare GmbH - Germany

Siemens Shenzhen Magnetic Resonance Ltd. - China

Siemens Shanghai Medical Equipment Ltd. - China

Siemens Healthcare Private Ltd. - India

Siemens Medical Solutions USA Inc. - USA

Siemens Healthcare Diagnostics Inc. - USA and/or

Siemens Healthcare Diagnostics Products GmbH - Germany

<b>1</b>	<b>Materials</b>	<b>5</b>
	1.1 Materials	5
<b>2</b>	<b>Customer Information</b>	<b>6</b>
	2.1 Site Survey Customer Information	6
	2.1.1 Institution	6
	2.1.2 Client Laboratory Contact	6
	2.1.3 Preliminary Visit	6
	2.1.4 Final Visit	6
	2.1.5 Site Surveyor	7
	2.1.6 Service Manager or Designee	7
<b>3</b>	<b>Delivery Requirements</b>	<b>8</b>
	3.1 Delivery Information	8
	3.1.1 Unloading Area	9
	3.1.2 Installation Route	10
<b>4</b>	<b>Environmental Requirements</b>	<b>12</b>
	4.1 Environmental Requirements	12
	4.1.1 System Specifications	12
	4.1.2 Site Requirements	13
	4.1.3 Customer Site	16
<b>5</b>	<b>Power Requirements</b>	<b>23</b>
	5.1 Power Requirements	23
	5.1.1 Site Requirements	23
	5.1.2 Customer Site	26
<b>6</b>	<b>Water Requirements</b>	<b>28</b>
	6.1 Water Requirements	28
	6.1.1 Site Requirements	28
	6.1.2 Customer Site	28
<b>7</b>	<b>Drain Requirements</b>	<b>30</b>
	7.1 Drain Requirements	30
	7.1.1 Customer Site	30
<b>8</b>	<b>Air Supply Requirements</b>	<b>31</b>
	8.1 Air Supply Requirements	31
	8.1.1 Air Supply Information	31
<b>9</b>	<b>Remote Connectivity Requirements</b>	<b>36</b>
	9.1 Remote Connectivity Requirements	36
	9.1.1 Customer Site	36

<b>10</b>	<b>Other Requirements</b>	<b>40</b>
10.1	Other Requirements	40
10.1.1	Actual Equipment Location	41
10.1.2	Sample Container Requirements	43
10.1.3	Automation Requirements	46
<b>11</b>	<b>Notes / Comments</b>	<b>48</b>
11.1	Notes / Comments	48
<b>12</b>	<b>Changes to Previous Version</b>	<b>49</b>
12.1	Changes to Previous Version	49
12.1.1	Version 13	49
12.1.2	Version 12	49
12.1.3	Version 11	50
12.1.4	Version 10	51
12.1.5	Version 09	51
12.1.6	Version 08	51
12.1.7	Version 07	51
12.1.8	Version 06	52
12.1.9	Version 05	52
12.1.10	Version 04	52
12.1.11	Version 03	53
12.1.12	Version 02	53

## 1.1 Materials

- Tape measure
- String and tape
- Laser level (if available)

## 2.1 Site Survey Customer Information

### 2.1.1 Institution

Address:			
City:		State/Province:	
Country:		Postal/Zone Code:	

### 2.1.2 Client Laboratory Contact

Name:			
Telephone:			
Date Reviewed Preliminary Report:			
Signature:			
Date Reviewed Final Report:			
Signature:			

### 2.1.3 Preliminary Visit

Date:			
Signature:			
Is Site Ready For Installation?	<input type="checkbox"/> Yes <input type="checkbox"/> No		

### 2.1.4 Final Visit

Date:			
Signature:			
Is Site Ready For Installation?	<input type="checkbox"/> Yes <input type="checkbox"/> No		

### 2.1.5 Site Surveyor

Name:	<input type="text"/>
Telephone:	<input type="text"/>

### 2.1.6 Service Manager or Designee

Signature:	<input type="text"/>	
Siemens Support Personnel		
<b>Name</b>	<b>Contact Information</b>	
<input type="text"/>	<input type="text"/>	
<input type="text"/>	<input type="text"/>	
<input type="text"/>	<input type="text"/>	
<input type="text"/>	<input type="text"/>	

### 3.1 Delivery Information



The information provided here is for planning delivery. For space requirements of an installed configuration in the laboratory, refer to (→ Space Requirements / Page 13).



During installation, ensure there is enough room in the laboratory to stage boxes from the delivery. These boxes contain system installation parts, including covers and parts that must be installed in an analyzer, DL, SH, SHC or Decapper.

Each configuration will have different numbers of boxes.

Tab. 1 System Dimensions\* (see note below)

	Height	Width	Depth	Weight
<b>Sample Handler</b>	1465 mm (57.7 in.)	876.3 mm (34.5 in.)	1431.3 mm (56.4 in.)	464.5 kg (1024 lb.)
<b>Chemistry Analyzer</b>	1363.5 mm (53.7 in.)	1452.5 mm (57.2 in.)	1183.3 mm (46.6 in.)	470.4 kg (1037 lb.)
<b>Immunoassay Analyzer</b>	1500 mm (59.1 in.)	1452.5 mm (57.2 in.)	1167 mm (45.9 in.)	594.7 kg (1311.1 lb.)
<b>Direct Load</b>	1365 mm (53.7 in.)	425.3 mm (16.7 in.)	1150 mm (45.3 in.)	124.3 kg (273 lb.)
<b>Sample Handler Connect</b>	1350 mm (53.1 in.)	400 mm (15.7 in.)	1350 mm (53.1 in.)	117.9 kg (260 lb.)
<b>Decapper</b>	1223 mm (48.1 in.)	254 mm (10 in.)	1109 mm (43.7 in.)	79.4 kg (175 lb.)



Dimensions given include:

- track for a track system (analyzers / SH / SHC / Decapper) but exclude the Magline loop track between the SH and an Analyzer
- reagent loader, panels and covers in place for a stand-alone system (analyzers / DL)

Tab. 2 Crated Analyzer / SH / SHC / DL / Decapper Dimensions

	Height	Width	Depth	Weight
<b>Sample Handler</b>	1752.6 mm (69.0 in.)	2235.2 mm (88.0 in.)	1270 mm (50.0 in.)	645.1 kg (1422 lb.)
<b>Chemistry Analyzer</b>				630 kg (1388.9 lb.)
<b>Immunoassay Analyzer</b>				770 kg (1697.6 lb.)
<b>Direct Load</b>	1651 mm (65.0 in.)	2260.6 mm (89.0 in.)	1041.4 mm (41.0 in.)	287.4 kg (633.6 lb.)
<b>Sample Handler Connect</b>				242.9 kg (535.5 lb.)
<b>Decapper</b>	1353 mm (53.25 in.)	572 mm (22.5 in.)	1181 mm (46.5 in.)	123 kg (271 lb.)

### 3.1.1 Unloading Area

1. Is there a loading dock?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
2. If a loading dock is present, is the loading dock height adjustable?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If the above answer is No, what is the height of the loading dock from the ground?		
3. If a loading dock is not present, special lifting equipment must be provided by the trucking company for a "sidewalk delivery." Will this be a sidewalk delivery?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If the above answer is Yes, the trucking company will require a description of the unloading area, including the distance from the unloading area to the building and any special lifting requirements. Provide a brief description of the unloading area in the Comments field below.		
Comments:		
4. Is there enough space on the loading dock or somewhere else to support the delivery of the crated module configuration (e.g., DL/CH or IM, Scl, Sccll, etc.) and perform the uncrating of the modules the site is ordering?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
5. Can a 14 m (46 ft.) trailer be backed into the loading area?	Yes <input type="checkbox"/>	No <input type="checkbox"/>



6. On which days and during what hours can the shipment be received?		
Comments:		
<input type="text"/>		
Preliminary Visit: Accepted <input type="checkbox"/>	Rejected <input type="checkbox"/>	
Comments:		
<input type="text"/>		
Final Visit: Accepted <input type="checkbox"/>	Rejected <input type="checkbox"/>	
Comments:		
<input type="text"/>		

### 3.1.2 Installation Route

Each analyzer / SH / SHC / DL / Decapper fits through a 1067 mm (42.0 in.) delivery pathway without the need to remove any subassemblies.

Each analyzer / SH fits through an 863 mm (34.0 in.) delivery pathway after the removal of some subassemblies.

1. Do all door turns on the installation route have a minimum radius of 1690 mm (66.5 in.)?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
Comments:		
<input type="text"/>		
2. Do all doors and hallways along the installation route meet the following requirements? <b>Minimum width:</b> 1067 mm (42.0 in)	Yes <input type="checkbox"/>	No <input type="checkbox"/>
Comments:		
<input type="text"/>		
3. Will an elevator be used during delivery?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If the above answer is Yes, please note the dimensions of the elevator (including the door opening dimensions) in the Comments field below.		
Comments:		
<input type="text"/>		
4. Are there stairs on the installation route?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If the above answer is Yes, please note the number of stairs and the width of the stairway(s) in the Comments field below.		
Comments:		
<input type="text"/>		
5. Is the floor along the delivery route smooth and free of sills/gaps?	Yes <input type="checkbox"/>	No <input type="checkbox"/>

If the above answer is No, ramps are required. Please indicate the number of ramps required in the Comments field below.		
Comments:		
6. Describe any other obstacles or factors that may impact delivery in the Comments field below.		
Comments:		
<b>Preliminary Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:		
<b>Final Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:		

## 4.1 Environmental Requirements

### 4.1.1 System Specifications

#### 4.1.1.1 Thermal Output Specifications

Tab. 3 Thermal Output

Type	Idle	Running
CH	1300 Watt 4450 BTU/hr	1530 Watt 5210 BTU/hr
DL	154 Watt 525 BTU/hr	300 Watt 1024 BTU/hr
IM	1300 Watt 4450 BTU/hr	1330 Watt 4530 BTU/hr
SH	735 Watt 2506 BTU/hr	1073 Watt 3661 BTU/hr
SHC	386 Watt 1317 BTU/hr	866 Watt 2955 BTU/hr
Decapper	240 Watt 819 BTU/hr	720 Watt 2456 BTU/hr

#### 4.1.1.2 Noise Specifications

The System shall meet applicable Noise Rating (NR) and Noise Criteria (NC) for medical devices in laboratories. Note: NR is used for European markets, NC is used for United States and Non-European markets.

The average sound pressure level for any analyzer / SH / DL / Decapper when utilizing the ISO3744 25 or 9 point measuring technique or the minimum of 9 point (1 meter away with analyzer / SH / DL / Decapper on a reflecting plane) shall satisfy NC-46. Exception: the IM analyzer shall satisfy NC-55.

The average sound pressure level for an Scl when utilizing the ISO3744 25 or 9 point measuring technique or the minimum of 9 point (1 meter away with analyzer / SH on a reflecting plane) shall satisfy NC-45.

Below is the average noise output from an SH, the analyzers, the Decapper and an Scl configuration:

- SH: 53 dBA
- CH: 50 dBA
- IM: 65 dBA
- Scl: 65 dBA

- Decapper: 46 dBA

## 4.1.2 Site Requirements

### 4.1.2.1 Space Requirements

The following clearances around an installed system are required (left and right are from the front perspective of the system):

- 609.6 mm (24 in.) on the right side of the system. (Maintain this clearance for installations that are being expanded, for example, with the addition of a Decapper at the end.)
- 609.6 mm (24 in.) on the left side of DL or non-connected SH configuration.
- 914 mm (36 in.) **minimum** / 1371.6 mm (54 in.) **recommended** in the front of a SH configuration.
- 609.6 mm (24 in.) in the front of a DL configuration.
- 914 mm (36 in.) in the rear of any system.

Refer to the following dimensions when planning the laboratory layout for the Atellica Solution configuration.



Adding all of the width values together for the applicable components / modules listed in (→ Tab. 1 Page 8) will most likely result in a width value larger than the actual configuration width. This is due to overlap of some components / modules when they are installed together. To calculate the dimension values of a specific configuration, use the EasyCalc tool (→ EasyCalc Tool / LDAT-000.840.21).



The overall width of a Decapper installed between analyzers or on the end of a configuration is the width of a Decapper plus 90 mm (3.5 in.) as a result of the overlap of the track on each side of the Decapper (approximately 45 mm (1.75 in.)). A Decapper in the Magline loop does not add any additional width to a configuration.

Below are the footprints of several common configurations for the Atellica Solution, including the applicable service areas.



The following images show the minimum service clearance.

Fig. 1: Atellica SH with Analyzers (Scl)




Fig. 2: Atellica Solution with Multiple Sample Handlers (SScII)




Fig. 3: Atellica SH with Analyzers and Sample Handler Connect (sScI)






Fig. 4: Atellica CH Analyzer with Direct Load



**i** The widths and depths are the same for an Atellica IM Analyzer with Direct Load as above.

Fig. 5: Atellica SH with Analyzers and Decapper (ScDI)



### 4.1.3 Customer Site

<p>1. The Atellica Solution System requires:</p> <ul style="list-style-type: none"> <li>■ Ventilation sufficient to maintain 18–30 °C (64–86 °F) operating temperature, with a maximum temperature change of 3 °C/hour.</li> <li>■ Relative humidity of 20–80% with system operating with no condensation.</li> </ul>	
<p>Has the assigned plant/maintenance engineer verified that throughout the year, the proposed system location meets these requirements?</p>	<p>Yes <input type="checkbox"/> No <input type="checkbox"/></p>

2. Space around the system must be sufficient to allow access by the user. There must be a minimum of:		
<ul style="list-style-type: none"> <li>■ 91.4 cm (36 in.) behind the system (connected, non-connected, and DL)</li> <li>■ 61.0 cm (24 in.) on the left of the system (non-connected system or DL system)</li> <li>■ 61.0 cm (24 in.) on the right of the system (connected, non-connected, and DL)</li> <li>■ 91.4 cm (36 in.) in the front of the system; 61.0 cm (24 in.) in front of a DL system.</li> </ul>		
Does the proposed system location meet these requirements?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
3. Air conditioning and heating ducts should not blow directly on the system, nor should the system be exposed to direct sunlight.		
Does the proposed system location take these restrictions into consideration?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
4. Did you discuss the potential for locating the system to allow space for future expansion with the customer?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
5. Is a Decapper being added between analyzers to an existing system installation? If the response is Yes, order new water and waste tubing to accommodate the additional length of the installation.	Yes <input type="checkbox"/>	No <input type="checkbox"/>
6. Is the site elevation between 0 m and 2000 m (6562 ft.)?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
What is the site's elevation in meters ( $\pm 100$ m)?		
7. The system should be placed on a level floor (45 mm maximum height variation) that is capable of supporting a load* specific for the customer's configuration that includes two or more of the following:	<b>*Floor Load</b>	
<b>SH:</b> 464.5 kg (1035 lb) / 1221.5 mm (48.1 in) W x 1431.3 mm (54.4 in) D =	266 kg/m <sup>2</sup>	
<b>DL:</b> 124.3 kg (273 lb) / 425.3 mm (16.7 in) W x 1150 mm (45.3 in) D =	254 kg/m <sup>2</sup>	
<b>CH:</b> 470.4 kg (962 lb) / 1452.5 mm (57.2 in) W x 1183.3 mm (46.6 in) D =	274 kg/m <sup>2</sup>	
<b>IM:</b> 594.7 kg (1308 lb) / 1452.5 mm (57.2 in) W x 1167 mm (46.9 in) D =	351 kg/m <sup>2</sup>	
<b>SHC:</b> 117.9 kg (260 lb) / 400 mm (15.7 in) W x 1350 mm (53.1 in) D =	218 kg/m <sup>2</sup>	
<b>Decapper:</b> 79.4 kg (175 lb) / 250 mm (9.8 in) W x 1109 mm (43.7 in) D =	286 kg/m <sup>2</sup>	
*The module floor load is presented above. There will be higher concentrated loads at each caster or leveling foot (see below).		
<b>Distance of Feet from Each Other**</b>	<b>Load per Foot</b>	



<b>SH:</b> 715 mm (side to side); 1050 mm (front to back)	122.5 kg
<b>DL:</b> 768.3 mm (front to back)	62.4 kg
<b>CH:</b> 1283.6 mm (side to side); 713.1 mm (front to back)	126.6. kg
<b>IM:</b> 1283.6 mm (side to side); 713.1 mm (front to back)	155.5 kg
<b>SHC:</b> 215.5 mm (side to side); 1050 mm (front to back)	29.5 kg
<b>Decapper:</b> 151.0 mm (side to side); 739.5 mm (front to back)	19.9 kg
**Refer to images in the "Distances Between Feet" section (→ Distances Between Feet / Page 19) for feet distances.	
Can the floor support the load of the configuration?	Yes <input type="checkbox"/> No <input type="checkbox"/>
8. Is seismic anchoring required?	Yes <input type="checkbox"/> No <input type="checkbox"/>
9. What is the pitch of the floor? Use the steps below this table to calculate the pitch and enter the value in the field below.	
<input type="text"/>	
<b>Preliminary Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>	
Comments:	
<input type="text"/>	
<b>Final Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>	
Comments:	
<input type="text"/>	

#### 4.1.3.1 Determine the Pitch of the Floor




This section is not necessary if using a laser level.

1. Attach a long piece of string (at least as long as the length of the Atellica Solution configuration) to a ruler.
2. Place the ruler-end of the string at the end of the configuration space that is suspected to be the lower of the 2 endpoints.
3. Place the other end of the string at the end of the configuration space that is suspected to be the higher of the 2 endpoints.
4. With a second person holding the end of the string on the floor, raise the string up on the ruler end until the string appears level.
5. With a third person measuring the string so that is level, continue to adjust the height of the string on the ruler until the string is level.
6. Measure the height of the string on the ruler.

If the height is greater than 1 1/4", 1 spacer per leveling foot (4 for each analyzer / SH / SHC / Decapper; 2 for DL) is required for each multiple of 1 1/4". For example, if the height is 3" on an Scl system, 8 spacers will be required for the IM analyzer.


#### 4.1.3.2 Distances Between Feet

Fig. 6: SH - Distances Between Feet




- (1) Side to Side: 715 mm
- (2) Front to Back: 1050 mm

Fig. 7: DL - Distances Between Feet



(1) Front to Back: 768.3 mm


Fig. 8: CH - Distances Between Feet



(1) Front to Back: 713.1 mm


(2) Side to Side: 1283.6 mm

Fig. 9: IM - Distances Between Feet




- (1) Front to Back: 713.1 mm
- (2) Side to Side: 1283.6 mm

Fig. 10: SHC - Distances Between Feet



- (1) Side to Side: 215.5 mm
- (2) Front to Back: 1050 mm

Fig. 11: Decapper - Distances Between Feet




- (1) Front to Back: 739.5 mm
- (2) Side to Side: 151.0 mm

## 5.1 Power Requirements

### 5.1.1 Site Requirements

- The Atellica Solution analyzers / SH / DL are shipped with a NEMA 6-20P (230VAC, 20A) plug. The plug can be changed if necessary to meet local requirements.
- Each analyzer or SH (not including SHC, DL, or DC, which is powered through the analyzer) requires a 4.4 kVA (USA) or 3.7 kVA (Europe) power source, single-phase, two-pole, three-wire configuration with Class III grounding. An analyzer / SH shall support incoming AC voltage from a nominal line voltage range of 200-240VAC, 50/60Hz. The main supply voltage fluctuations are not to exceed  $\pm 10$  percent of the nominal voltage.
- Each Decapper requires a .72 kVA (USA) or .69 kVA (Europe) power source, single-phase, two-pole, three-wire configuration with Class III grounding. A Decapper shall support incoming AC voltage from a nominal line voltage of 240v, 50/60Hz. The main supply voltage fluctuations are not to exceed  $\pm 10$  percent of the nominal voltage.
- The Atellica Decapper is shipped with an IEC 60320 C14 recessed power inlet. The Decapper requires a 5.0 meter power cord with an IEC 60320 C13 connector and a plug suitable for use by local regulations that matches the installed receptacle.

Tab. 4 Decapper Power

Con- nector	Appli- ance In- let	Configuration (Connector/Inlet)	Groun- ded	Max Current		Voltage	
				Int'l	N.Amer- ica	Int'l	N.Amer- ica
C13	C14	<p><i>Fig. 12: C13 Plug Connector/C14 Power Inlet</i></p> 	Yes	10A	15A	250V	125/250 V


- The Atellica Solution system conforms to Installation/Overvoltage Category II (as specified in IEC 61010-1).
- The Maximum Inrush Current is 250 A @ 3.3 ms for Atellica Solution analyzers.
- The circuit must have a separate, dedicated (NOT ISOLATED) line L1, L2 and dedicated (NOT ISOLATED) ground in its own conduit. The conduit should start at the distribution panel and be continuous to the receptacle. Three-wire distribution to the receptacle is required for each Atellica Solution analyzer, DL, SH, or SHC. The third (green or green/yellow) ground wire should start at the distribution panel and be continuous to the receptacle in accordance with NEC paragraph 250.146 (D), exception 4, unless local codes prohibit. The ground wire can not be tied to grounds from other loads.
- The customer must provide a receptacle installed by a qualified electrician before arrival of the Atellica Solution. The receptacle must be accessible to the 5.0 meter power cord furnished with the Atellica Solution. The USA requires 20 Ampere, 250 Volt receptacle, Hubbell receptacle # HBL5461. Outside USA requires a 16 Ampere or greater receptacle.
- Where the specified North America plug cannot be used, the customer must provide a power plug suitable for use by local regulations that matches the installed receptacle.

The power plug must be compatible with 3 conductor 3.31 mm<sup>2</sup> and the strain relief able to accept minimum 11.2 ± 0.15 mm diameter cordage.


Tab. 5 Current Draw

Mode	Current (60 Hz, 220 V)	Current (50 Hz, 220 V)
In-Rush with Main Breaker ON	250 A @ 3.3 ms	250 A @ 3.3 ms
Operation- Peak	8 A	8 A
Main Circuit Breaker Rating	20 A long delay in USA, outside USA use MCB 16 A Type-C or 16 A Type-D	Outside USA use MCB 16 A Type-C or 16 A Type-D.

Tab. 6 Mains Input Plug

Cat. No. (Black Cover) (Black Cord Clamp)	Cat. No. (Clear Cover) (Natural Cord Clamp)	Cat No. Hospital Grade (Clear Cover) (Nat. Cord Clamp)	Configuration	Rating	UL	CSA
HBL5466CA			Fig. 13: NEMA 6-20P 	20A / 250V	Yes	Yes

Tab. 7 Mains Receptacle

Catalog Numbers			Hospital Grade	Color	Configuration	Rating
Ring Terminal	Standard	Marine				
	HBL5461			Brown	Fig. 14: NEMA 6-20R 	20A / 250V
	HBL5161I			Ivory		

#### 5.1.1.1 External UPS (Optional) Requirements

New system installations that will use UPS should be wired using a circuit and receptacle to match the specifications of the UPS. If the system will not be using a UPS, then a 20 A (USA) / 16 A (Outside USA) circuit per analyzer / SH / DL / SHC is sufficient.

Tab. 8 External UPS Minimum Requirements: Input / Output

Input	Value	Output	Value
Voltage (VAC)	200 to 240	Capacity	3000VA / 2700W
Phase	Single	Voltage (VAC)	200 / 208 / 220 / 230 / 240
Frequency (Hz)	50 / 60	Voltage Regulation	±3.0% Max

Input	Value	Output	Value
Input Connection	NEMA L6-20P	Output Voltage THD	< 3.0% for linear loads < 5.0 % for non-linear loads
Generator Compatibility	Yes	Frequency (Hz)	50 / 60
		Frequency Regulation	± 0.5 Hz
		Output Connections	L6-30R L6-20R

Tab. 9 External UPS Minimum Requirements: Battery

Attribute	Value
Back Up Time (minutes)	≥ 8.0 @ Full Capacity
Recharge Time	≤ Recharge to 90% within 10x's actual discharge time

Tab. 10 External UPS Minimum Requirements: Environment

Attribute	Value
Operating Temperature (°C)	10 to 40
Operating Altitude (meters)	2000
Humidity (% RH non-condensing)	5 to 90
Heat Dissipation (BTU/Hour)	< 1650

Tab. 11 External UPS Minimum Requirements: Agencies

Attribute	Value
EMC	CISPR Class B FCC part 15 Class B / ICES-003 IEC61000-4-2 IEC61000-4-3 IEC61000-4-4 IEC61000-4-5 IEC61000-4-6 IEC61000-4-8
Safety	IEC 62040-1:2008



Tab. 12 External UPS Minimum Requirements: Other

Attribute	Value
Topology	True On-Line, Double Conversion, Pure Sine Wave
Communication	None
Seismic Compatibility	Yes

### 5.1.1.2 Compressor (Decapper-specific) Requirements

Installations that will use a Decapper and do not have house air will need a compressor. The installation should be wired using a circuit and receptacle to match the specifications in the compressor documentation.

### 5.1.2 Customer Site

1. List the system configuration(s) to be installed and the total number of dedicated circuits required (Examples: SCI – 3, DL+ IA – 1):		
2. Will this installation include an external UPS for each analyzer or SH?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If Yes, please list the total number of UPSs and for which analyzers / SH / DL:		
3. Do the AC power outlets meet the specifications listed in the Power Requirements section above?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If the AC power outlets do not meet the specifications above, instruct the customer to meet with the facility engineer(s) to ensure AC power outlets are installed that can meet the specifications.		
4. Siemens Healthineers recommends a dedicated circuit for each analyzer or SH with or without a UPS. Are dedicated circuits available that meet the circuit breaker ratings below?		
<b>US:</b> 220 V, 20 A <b>Europe:</b> 230 V, 16 A Type C or Type D	Yes <input type="checkbox"/>	No <input type="checkbox"/>
5. Siemens Healthineers recommends a dedicated power receptacle for each Decapper. Is there a dedicated power receptacle available that meets the ratings below?		
<b>US:</b> 220 V, 3 A <b>Europe:</b> 230 V, 3 A	Yes <input type="checkbox"/>	No <input type="checkbox"/>
6. Will a compressor be needed to operate the Decapper?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If the above answer is Yes, the customer must supply a power receptacle for the compressor.		

7. List the number of power receptacles available for the Atellica system.		
8. Is an extra power receptacle for a printer located near the system and available?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
9. The power cables are 5 m (16 ft.) long. Will these cables be routed across walkways?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If the above answer is Yes, the customer must supply cable guards, as well as a power receptacle.		
10. If an emergency generator is used, is a power conditioner available?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
11. In the space below, provide a diagram of the AC receptacle, or list the AC power outlet configuration.		
Comments:		
<b>Preliminary Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:		
<b>Final Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:		

## 6.1 Water Requirements

### 6.1.1 Site Requirements

- Special Reagent Grade water (SRW) (CLSI C3-A4 Reagent Laboratory Water, with exceptions noted below).
- Individual water supply units can be used or a single, centrally located system can supply the water to all analyzers.
- IM and CH analyzers accept incoming reagent grade water from the customer facility at a regulated pressure between 5 PSI and 30 PSI.
- The CH analyzer consumes less than 33 liters of water per hour for steady state operation.
- The IM analyzer consumes less than 6 liters of water per hour for steady state operation.
- The IM and CH analyzers demand incoming water at a flow rate between 2.0 L/min. and 4.5 L/min (depending on the pressure).
- Temperature of water from the supply: 10-30 °C
- The system is supplied with 10 m (33 ft.) of 3/8" OD x .245 ID Polyurethane tubing, which is used for the SRW water connection. Size the SRW water fittings accordingly and use only compression type fittings due to water hammer effects.
- Acceptable CFU/mL (by test source):
  - Laboratory water system output  $\leq$  50 CFU/mL\*
  - \* Based on CLSI/CAP Reagent Grade water specifications

### 6.1.2 Customer Site

1. Is SRW located on site?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If the above answer is No, a SRW supply system able to supply pressurized water will be required.		
Comments:		
2. Does the water meet the following specifications for Special Reagent Water (SRW)?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<ul style="list-style-type: none"> <li>■ Resistivity (<math>\geq</math> 10 M<math>\Omega</math>-cm)</li> <li>■ Bacteria (<math>\leq</math> 50 CFU/mL)</li> <li>■ TOC – Total Organic Carbon (<math>\leq</math> 500 ppb)</li> <li>■ A 0.22 micron filter is required at the output stage of the laboratory water purification system and an additional 0.22 micron filter is required before the input to the water supply for each analyzer</li> </ul>		
Refer to the Water Quality Testing document ( $\rightarrow$ Water Quality Testing / LDAT-000.840.10) for testing information.		

Comments:		
3. Some laboratories may require special pre-filters to the SRW unit. Are filters required? For example:	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<ul style="list-style-type: none"> <li>▪ Sand and carbon filters for dissolved organics</li> <li>▪ Colloidal filter for high iron content water</li> <li>▪ Pretreatment (softening) kit for dissolved inorganics</li> </ul> Consult with a local water purification company if necessary. List any required filters in the Comments field below.		
Comments:		
4. Are the SRW parameters within the specifications listed in question 2 above? List any parameters that do not meet the specifications in the Comments field below.	Yes <input type="checkbox"/>	No <input type="checkbox"/>
Comments:		
5. Is the water pressure regulated to between 5 PSI (min) and 30 PSI (max)?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
Comments:		
<b>Preliminary Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:		
<b>Final Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:		

## 7.1 Drain Requirements

### 7.1.1 Customer Site

1. Is a floor drain located within 9.1 m (30 ft.) of the back of the system?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
2. A drain can be located at a maximum height of 122 cm (48 in) from the floor. Is the drain height acceptable?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
3. Is the drain capable of handling a minimum of 100 L/hr. (1.7 L/min) per analyzer?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
4. The waste pipe inlet must be large enough to accommodate various tubes. A 76.2 mm (3 in.) diameter drain pipe is recommended.		
Is the diameter of the drain pipe acceptable?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
5. Do the local laboratory practices and/or applicable environmental regulations allow for waste to be routed into the laboratory's drain?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
6. If a compressor requiring a drain is installed to support a Decapper, is the drain located within the distance specified in the compressor documentation?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If the above answer is No, are there plans to install a drain?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If the response to the previous question is No, the drainage system will require an upgrade.		
7. On what date will this upgrade be completed?	Date <input type="text"/>	
<b>Preliminary Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:	<input type="text"/>	
<b>Final Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:	<input type="text"/>	

## 8.1 Air Supply Requirements

The Decapper has no internal air compressor and requires the use of either laboratory air supply or a client-supplied air compressor.

### 8.1.1 Air Supply Information

#### 8.1.1.1 System Requirements

The Atellica Solution Decapper module:

- consumes 1.76 NLPM (.062 CFM) of compressed air per minute
- provides air cut-off and air regulator control with an air pressure gauge on the rear panel
- has a 6 mm pneumatic tubing input connection on the rear panel
- should, once connected to the air supply, show an air pressure reading between .6 MPa (Megapascal) and .65 MPa (87 PSI to 94 PSI) on the gauge at the back of the Decapper
- operates at an altitude of up to 2000 meters (6562 ft.)

#### 8.1.1.2 Site Requirements

Either laboratory air supply or client-supplied air compressor that:

- supplies between 90 and 100 PSI to a Decapper
- supplies an air flow of 1.76 NLPM (.062 CFM) per minute to a Decapper
- uses a compressed air dryer
- uses a drain for condensation run-off
- meets the following compressed Air Purity Class 1 output rating (clean, oil-free and dry compressed air) per ISO 8573-1:2010:

Class Item	Values
Particles	<ul style="list-style-type: none"> <li>■ Less than or equal to 20,000 particles per cubic meter (between 0.1 micrometer and 0.5 micrometer)</li> <li>■ Less than or equal to 400 particles per cubic meter (between 0.5 micrometer and 1.0 micrometer)</li> <li>■ Less than or equal to 10 particles per cubic meter (between 1.0 micrometer and 5.0 micrometer)</li> </ul>
Humidity and Liquid Water	Pressure dew point less than or equal to -70 degrees Celsius
Total Oil	Less than or equal to 0.01 milligrams per cubic meter

## 8.1.1.3 Customer Site

1. Does the laboratory have a single, centrally located air supply system? If the response is No, proceed to question 10.	Yes <input type="checkbox"/>	No <input type="checkbox"/>
2. Does the system have a pressure regulator, pressure gauge and control valve to ensure appropriate air supply to one or more Decappers? If the response is No, these items should be installed.	Yes <input type="checkbox"/>	No <input type="checkbox"/>
3. Does the air supply have an air filtration system that meets the specified Air Purity Class 1 output rating? If the response is No, a filter/dryer needs to be installed.	Yes <input type="checkbox"/>	No <input type="checkbox"/>
4. Each Decapper consumes 1.76 NLPM (.062 CFM) per minute. Enter the total required rate of air flow, in either Normal Liters per minute (NLPM) or Cubic Feet per minute (CFM).	NLPM <input type="text"/> CFM <input type="text"/>	
5. What is the current CFM usage of the air supply system? What is the minimum CFM supplied by the air supply system?	Current CFM <input type="text"/> Minimum CFM <input type="text"/>	
6. The required air supply to be delivered to each Decapper is between .6 MPa (Megapascal) and .7 MPa (87 PSI to 101.5 PSI). Confirm that the required pressure is available at each Decapper location.	Yes <input type="checkbox"/>	No <input type="checkbox"/>
7. What is the minimum pressure provided by the air supply system?	Min PSI <input type="text"/>	
8. Can the existing air supply system handle the additional air requirements (CFM and PSI) listed above?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
Comments: <input type="text"/>		
If No, are there plans to improve the air supply system?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If the response to the previous question is No, the house-supplied air supply system will require an upgrade.		
9. On what date will this upgrade be completed? Go to question 21.	Date <input type="text"/>	
10. <b>Atellica Systems requiring an air compressor unit:</b> Each Decapper consumes 1.76 NLPM (.062 CFM) per minute. Enter the total required rate of air flow, in either Normal Liters per minute (NLPM) or Cubic Feet per minute (CFM).	NLPM <input type="text"/> CFM <input type="text"/>	

<p>11. The required air supply to be delivered to each Decapper is between .6 MPa (Megapascal) and .7 MPa (87 PSI to 101.5 PSI). Confirm that the required pressure is available at each Decapper location.</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>12. Will the air compressor be operating at a high altitude? If Yes, compressor performance at higher altitudes needs to be taken into account.</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>13. Does the air supply have an air filtration system that meets the specified Air Purity Class 1 output rating? If the response is No, a filter/dryer needs to be installed.</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>14. Did you provide the laboratory supervisor with CFM, PSI and altitude requirements for compressor sizing?</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>15. Did you discuss ordering requirements and delivery timing with the laboratory supervisor?</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>16. Establish a date on which to follow up with the laboratory supervisor to confirm the compressor order.</p>	Date _____	
<p>When will the compressor be on site?</p>	Date _____	
<p>17. Did the laboratory supervisor provide a date on which the location designated for the compressor will be evaluated to confirm the installation requirements are met? Confirm that the location will provide sufficient air flow and pressure to the Decapper for proper functioning. Refer to the compressor technical specifications for additional information. NOTE: A locally licensed electrician will be needed to make AC power connections.</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>18. Are a gauge and shut-off valve with standard connecting thread installed?</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>19. Does the compressor have an internal accumulator/air receiver tank? If the response is No, does the customer have space for the tank? Confirm space for the tank is accounted for and that servicing of the tank is handled by the facility.</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>20. Is the plumbed line to specification? Refer to the compressor technical specifications for additional information.</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>21. Are any manifolds being installed to divide the air lines for the number of Decappers being installed? If Yes, factor in the placement of the manifold(s) into the answers for the path, air line length and cut off value questions.</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>22. What is the path of the air line between the air supply and</p>		



Decapper 1?		
Decapper 2?		
Decapper 3?		
23. What is the length of the air line needed for		
Decapper 1?	Length	
Decapper 2?	Length	
Decapper 3?	Length	
NOTE: It is good practice to install tubing a little longer than needed as the end could be cut off to ensure proper fit.		
24. Each Decapper needs a pneumatic cut off valve installed on the line. How many pneumatic cut off valves are needed?	Qty	
When will the valve be installed?	Date	
25. Have the cut off valves and the air lines been ordered?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
26. Have the cutoff valves and the air lines arrived on site?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
Comments:		
<b>Preliminary Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:		
<b>Final Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		

Comments:


## 9.1 Remote Connectivity Requirements

### 9.1.1 Customer Site



Serial connection or transmission control protocol/Internet protocol (TCP/IP) to Serial conversion is not supported. A TCP/IP network connection directly to the Laboratory Information System (LIS) Host must be used.

1. Is there an Atellica CM / sLCM already installed or planned for installation?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If No, then refer to the LCM CB-DOC DCIN-F00.811.06.XX.02 Pre-Installation and perform the Pre-Installation for the LCM while performing the Site Survey for the system.		
If Yes, record the Serial Number and IP Addresses (site and instrument) of the Atellica CM / sLCM.		
<b>Serial number:</b>		
<input type="text"/>		
<b>Site IP Address:</b>		
<input type="text"/>		
<b>Instrument IP Address:</b>		
<input type="text"/>		
If Yes, record the software version of the ACM / sLCM. <b>Note:</b> VA20B is the minimum software version that will support Atellica.		
<input type="text"/>		
2. What is the length of Network cable (CAT5) required to connect the SH Prime network port (right side rear of SH) to the hospital (WAN) network wall jack (e.g. 7ft., 10ft., 100ft.)?		
<input type="checkbox"/> 7-foot	10388887 CABLE, CAT V NETWORK 7 FOOT	
<input type="checkbox"/> 10-foot	10388888 CABLE, CAT V NETWORK 10 FOOT	
<input type="checkbox"/> 25-foot	10388889 CABLE, CAT V NETWORK 25 FOOT	
<input type="checkbox"/> 50-foot	10287529 CABLE, CAT V NETWORK 50 FOOT	
<input type="checkbox"/> 100-foot	10388890 CABLE, CAT V NETWORK 100 FOOT	
Ensure you order the appropriate length cable prior to installation of the Atellica Solution (not supplied with the system).		
3. Are there any existing middleware products in use or planned to be used with Atellica? If so, please complete the fields below with required information, or "N/A" if not applicable.		
<b>Middleware Products</b> (e.g Centralink, Inpeco DMS, OneLink, OEM Vendor):		
<input type="text"/>		

Local or Regional Automation/IT Consultant (e.g. Siemens or OEM Support Person):	
Local or Regional Automation/IT Consultant Email Address and Phone Number:	
Local or Regional Automation/IT Consultant Signature and Date:	
Signature:	Date:
4. Provide the following information about the site LIS:	
LIS Model	
IP Address	
Port	
Subnet Mask	
Default Gateway	

System Network Requirement/Permission	Responsible Person	Project Completion Date
Siemens requires a <b>dedicated network port jack</b> for the Atellica Solution to access Atellica CM / sLCM and Smart Remote Services (SRS). <b>Note: The Atellica Solution network uses a 192.168.1.X schema. If the customer's network uses this schema, the customer's IT department needs to change the schema so that the Atellica Solution can connect to the site's network.</b>		
Confirm the data being monitored by Siemens from the workstation with both the hospital lab manager and the hospital IT group and ensure that permission to monitor instrument data has been approved.		
Laboratory Manager:	Phone Number:	
Hospital IT/Security:	Phone Number:	
Comments:		

Site Network Connection Details	Responsible Person	Project Completion Date
<p><b>NOTE:</b> Static IP address assignment is preferred. If the customer prefers this method of connectivity, complete the fields below.</p> <p>Note: if static, it can not be a 192.168.1.x schema.</p>		
Complete if Static IP address will be provided.		
Yes <input type="checkbox"/> No <input type="checkbox"/>		
IP Address:	Preferred DNS Server:	
Subnet Mask:	Alternate DNS Server:	
Default Gateway:		
Will we require a proxy to access the Siemens Data Center?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If "Yes," record proxy server address:	Address:	Port:
Does the proxy require authentication?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If "Yes," add Username.	Username:	Password:
Are you using NTLM authentication (when using an ISA server as a proxy server)?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
Ensure that a Network port / jack (RJ45) connection will be available at the chosen location for the Atellica Solution. (Reserve and label the network jack for Siemens' use.)		
The internal network router for the Atellica Solution requires an administrator login and password. Record the customer's desired login and password below:		
Router Administrator Login:		
Router Administrator Password:		
Will the Atellica CM / sLCM support connectivity to the Atellica Solution over the customer network?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<b>If the above answer is Yes, the customer must allow access on common IPSEC ports (4500,500) to accommodate the secure connection</b>		

Site Network Connection Details	Responsible Person	Project Completion Date
Comments:		
<input type="text"/>		
<b>Preliminary Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:		
<input type="text"/>		
<b>Final Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:		
<input type="text"/>		

## 10.1 Other Requirements

<p>1. Inform the laboratory supervisor that the Atellica Solution System has the following barcode reading capabilities:</p> <ul style="list-style-type: none"> <li>▪ Maximum of 20 characters</li> <li>▪ I2 of 5*, Codabar*, Code 39*, Code 128</li> </ul> <p>* If enabled, I2 of 5 can use check digits and be sent to the host.</p>		
<p>Did you discuss the barcode capabilities with the laboratory supervisor?</p> <p>List the barcode symbologies to be used, as well as any additional required functionality:</p> <p>_____</p> <p>_____</p> <p>_____</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>2. Inform the laboratory about sample containers available for use with the Atellica Solution System.</p>		
<p>Did you discuss the available sample containers for use with the Atellica Solution System?</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>3. Does the customer provide support for laboratory staff to have personal mobile device access to their network using a wireless access point?</p> <p><b>Note:</b> The Atellica mobile web service requires the laboratory to provide their own wireless network access to use the feature.</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>4. Is cell phone usage permitted in the lab where the Atellica Solution is being installed?</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>If yes, is there cell phone coverage in the lab?</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>5. Is there an external printer already installed on the hospital network (WAN) to print reports for the Atellica Solution?</p>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
<p>If Yes, record the printer model type and number (e.g. Lexmark MS310DN-110) and network information.</p>		
<p>Printer Model</p> <p>_____</p>		
<p>IP Address</p> <p>_____</p>		
<p>Subnet Mask</p> <p>_____</p>		
<p>Default Gateway</p> <p>_____</p>		


6. Is there an Barcode label printer already installed on the hospital network (WAN) to print barcode labels for Atellica?	Yes <input type="checkbox"/>	No <input type="checkbox"/>
If Yes, record the printer model type and number (e.g. SATO CG212TT) and network information.		
Printer Model		
IP Address		
Subnet Mask		
Default Gateway		
<b>Preliminary Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:		
<b>Final Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:		

### 10.1.1 Actual Equipment Location

Refer to (→ Space Requirements / Page 13) for the recommended installation layout and space requirements. Use the space below to sketch a proposed location of the system in the laboratory. Indicate the approximate room dimensions, doorways and door swing, and distances from the system to the water supply, liquid waste drain and electrical supply.



Fig. 15: Actual Equipment Location



## 10.1.2 Sample Container Requirements

Tab. 13 Sample Container Requirements - Primary Tubes

Item	Material	Cap Type	Size	Volume	Note if false bottom shape	Tube Support by Configuration <sup>1</sup>		
						SH	DL	SHC
Plain Top Tubes ranging 13-16 dia, 75-100 high	plastic or glass	None	Various	4-10 mL	NA	Y	Some	N
13mm x 75 mm primary tube	plastic or glass	None	13mm x 75mm	4-5 mL	NA	Y	Y	Y
13x75 Vacutainer® Hemogard - BD	plastic	safety	13mm x 75mm	3.5 mL	NA	Y	Y	Y
13x75 Vacutainer - BD	plastic	stopper	13mm x 75mm	3.5 mL	NA	Y	Y	Y
13x75 Vacuette® - Greiner	plastic	safety	13mm x 75mm	4 mL	NA	Y	Y	Y
13mm x 100mm primary tube	plastic or glass	None	13mm x 100mm	6-7 mL	NA	Y	Y	Y
13x100 Vacutainer Hemogard - BD	plastic	safety	13mm x 100mm	6 mL	NA	Y	Y	Y
13x100 Vacutainer - BD	plastic	stopper	13mm x 100mm	6 mL	NA	Y	Y	Y
13x100 Vacuette - Greiner	plastic	safety	13mm x 100mm	5 mL	NA	Y	Y	Y
16mm x 100mm primary tube	plastic or glass	None	16mm x 100mm	10 mL	NA	Y	Y	Y
16x100 Vacutainer Hemogard - BD	plastic	safety	16mm x 100mm	8.5 mL	NA	Y	Y	Y
16x100 Vacutainer glass - BD	glass	stopper	16mm x 100mm	10 mL	NA	Y	Y	Y
16x100 Vacutainer - BD	plastic	stopper	16mm x 100mm	10 mL	NA	Y	Y	Y
16x100 Vacuette - Greiner	plastic	safety	16mm x 100mm	9 mL	NA	Y	Y	Y

Item	Material	Cap Type	Size	Volume	Note if false bottom shape	Tube Support by Configuration <sup>1</sup>		
						SH	DL	SHC
Sarstedt S-Monovette 16mm x 92mm	plastic	screw cap	16mm x 92mm	9 mL	NA	Y	N	Y
Sarstedt S-Monovette 13mm x 90mm	plastic	screw cap	13mm x 90mm	4.9 mL	NA	Y	Y	Y
Sarstedt S-Monovette 15.3mm x 75mm	plastic	screw cap	15.3mm x 75mm	5.5 mL	NA	Y	N	Y
Sarstedt S-Monovette 15mm x 92mm	plastic	screw cap	15mm x 92mm	7.5 mL	NA	Y	Y	Y
Trinidad Calibration / QC Tube 12x75 glass tube	glass-borosilicate	none	12mm x 75mm	5 mL	NA	Y	Y	Y
Bio-rad Control Tubes	plastic	thread	NA	4 mL	NA	Y	Y	N

1. The Decapper supports plastic tubes with safety and screw (thread) cap types and round and syringe bottom shapes. Tubes not supported by the Decapper can be used on the system uncapped. The operator's guide covers specific tube support.

Tab. 14 Sample Container Requirements - Transfer Tubes

Item	Material	Cap Type	Size	Volume	Bottom Shape	Tube Support by Configuration		
						SH	DL	SHC
Small Transfer Tube (12mm)	plastic or glass	None	12mm x 75mm	5 mL	Round	Y	Y	Y
Large Transfer Tube (13mm)	plastic or glass	None	13mm x 100mm	9 mL	Round	Y	Y	Y
Large Transfer Tube (14mm)	plastic	None	14mm x 100mm	5 mL	Round	Y	N	Y
Large Transfer Tube (16mm)	plastic or glass	None	16mm x 100mm	10 mL	Round	Y	Y	Y

						Tube Support by Configuration		
Item	Material	Cap Type	Size	Volume	Bottom Shape	SH	DL	SHC
Large Transfer Tube (16mm)	plastic	None	16mm x 75mm	8 mL	Round	Y	N	Y
Transfer Tube (Sarstedt) <sup>1</sup>	plastic	screw cap	15.3mm x 92mm	10 mL	Round	Y	N	Y

1. Supported by Decapper

Tab. 15 Sample Container Requirements - Sample Cups

						Tube Support by Configuration		
Item	Material	Cap Type	Size	Volume	Bottom Shape	SH	DL	SHC
Atellica TTSC - Siemens 11069061	plastic	None	13mm x 26mm	1 mL	Conical	Y	Y	N
Atellica TTSC - Siemens 11069062	plastic	None	13mm x 30mm	2 mL	Conical	Y	Y	N

Vista Tubes not on Atellica: BD Vacutainer® Conventional Stopper (10mm x 65mm and 10mm x 50mm)

Tab. 16 Capillary Collection Containers

					Tube Support by Configuration		
Item	Material	Cap Type	Volume	Bottom Shape	SH	DL	SHC
Becton Dickinson Microtainer® 365967	plastic	stopper	1 mL	Conical	Y	N	N
Greiner MiniCollect® 450472	plastic	stopper	1 mL	Conical	Y	N	N
Greiner MiniCollect Complete	plastic	stopper	0.8 mL	Conical (false-bottom)	Y	N	N
Sarstedt Microvette® 20.1344	plastic	stopper	1 mL	Conical	Y	N	N

Tab. 17 False Bottom Tubes

Item	Material	Cap Type	Size	Volume	Bottom Shape	Tube Support by Configuration		
						SH	DL	SHC
15.3 x 92 False Bottom - Sarstedt 62.612	plastic	screw	15.3mm x 92mm	4 mL	Flat	Y	Y	N

### 10.1.3 Automation Requirements

Tab. 18 Automation Software/Firmware Versions Needed for Atellica Sample Handler Connect Connectivity to Aptio/FlexLab

Automation System		Minimum Software/Firmware Version
Aptio by Siemens (Aptio-1)		Aptio 9.0
Aptio by Inpeco (Aptio-2)	SMS (automation track)	Quebec 1 update
	DMS (Inpeco Middleware)	Jena 1.0.1 update
FlexLab	SMS (automation track)	Quebec 1 update
	DMS (Inpeco Middleware)	Jena 1.0.1

1. Does the site meet the above automation requirements for use with Atellica Sample Handler Connect? <sup>1</sup>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
Comments:		
2. Will the Atellica system be connected to Aptio 1 or Aptio 2?	Aptio 1 <input type="checkbox"/>	Aptio 2 <input type="checkbox"/>
<b>Preliminary Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		
Comments:		
<b>Final Visit:</b> Accepted <input type="checkbox"/> Rejected <input type="checkbox"/>		

Comments:

- 
1. Atellica Decapper support is not provided with automation systems.

### 11.1 Notes / Comments



Multiple horizontal light gray bars representing a text input area for notes or comments.

## 12.1 Changes to Previous Version

### 12.1.1 Version 13

Change No.	Changes
LFT0012830	<b>4. Environmental Requirements</b> <ul style="list-style-type: none"> <li>(→ Customer Site / Page 16): Added a question on whether a Decapper is being added between analyzers to an existing installation, and if the response is Yes, to order new water and waste tubing.</li> </ul>
LFT0012929	<b>5. Power Requirements</b> <ul style="list-style-type: none"> <li>(→ Power Requirements / Page 23): Corrected the amperage value for Europe for the Decapper dedicated power receptacle.</li> </ul>
LFT0012703	<b>8. Air Supply Requirements</b> <ul style="list-style-type: none"> <li>(→ Air Supply Requirements / Page 31): Added a question on whether an in-house air supply system has a pressure regulator, pressure gauge and control valve to ensure appropriate air supply to one or more Decappers.</li> </ul>
LFT0013183	<b>10. Other Requirements</b> <ul style="list-style-type: none"> <li>(→ Other Requirements / Page 40): Added a statement to list the barcode symbologies to be used, as well as any additional required functionality, as part of the discussion with the laboratory supervisor about barcode capabilities.</li> </ul>
LFT0012940	<b>Throughout</b> <ul style="list-style-type: none"> <li>Corrected check mark functionality.</li> <li>Corrected broken links.</li> </ul> <p>No technical changes were made as part of these corrections.</p>

### 12.1.2 Version 12



This document supersedes LDAT-000.835.02.01.02.

Change No.	Changes
LFT0011548	<b>3. Delivery Requirements</b> <ul style="list-style-type: none"> <li>(→ Tab. 2 Page 9): updated all crated system weights except Decapper.</li> </ul>
LFT0011695 CHG0046330	<b>4. Environmental Requirements</b> <ul style="list-style-type: none"> <li>(→ Space Requirements / Page 13): added a note to this section regarding the EasyCalc tool, removed the configuration table because of the EasyCalc tool, and added updated and additional images for various Atellica Solution configurations.</li> </ul>



Change No.	Changes
LFT0011970	<p><b>5. Power Requirements</b></p> <ul style="list-style-type: none"> <li>(→ Power Requirements / Page 23): Added compressor power information.</li> </ul> <p><b>7. Drain Requirements</b></p> <ul style="list-style-type: none"> <li>(→ Drain Requirements / Page 30): Added compressor drain information.</li> </ul> <p><b>8. Air Supply Requirements</b></p> <ul style="list-style-type: none"> <li>(→ Air Supply Requirements / Page 31): Updated the value of the air pressure reading required for the Decapper.</li> </ul>
LFT0011326	<p><b>9. Remote Connectivity Requirements</b></p> <ul style="list-style-type: none"> <li>(→ Remote Connectivity Requirements / Page 36): added a note at the top of the section.</li> </ul>

### 12.1.3 Version 11

Change No.	Changes
CHG0050758	<p><b>3. Delivery Information</b></p> <ul style="list-style-type: none"> <li>(→ Delivery Information / Page 8): Updated to include Decapper specific delivery requirements.</li> </ul> <p><b>4. Environmental Requirements</b></p> <ul style="list-style-type: none"> <li>System Specifications (→ System Specifications / Page 12): Updated to include Decapper specific system requirements.</li> <li>Site Requirements (→ Site Requirements / Page 13): Updated to include Decapper specific site requirements</li> <li>Customer Site (→ Customer Site / Page 16): Updated to include Decapper specific customer site requirements.</li> </ul> <p><b>5. Power Requirements</b></p> <ul style="list-style-type: none"> <li>(→ Power Requirements / Page 23): Updated to include Decapper specific power requirements.</li> </ul> <p><b>8. Air Supply Requirements</b></p> <ul style="list-style-type: none"> <li>(→ Air Supply Requirements / Page 31): Added air supply requirements section for Decapper installations.</li> </ul> <p><b>10. Other Requirements</b></p> <ul style="list-style-type: none"> <li>Sample Container Requirements (→ Sample Container Requirements / Page 43): Added general information on Decapper tube specifications.</li> </ul>

### 12.1.4 Version 10

Change No.	Changes
LFT0010052	<b>8. Remote Connectivity Requirements</b> <ul style="list-style-type: none"> <li>(→ Remote Connectivity Requirements / Page 36): Updated the section to add notes related to potential conflicts between the customer and Atellica Solution network schemas.</li> </ul>
CHG0055019	<b>9. Other Requirements</b> <ul style="list-style-type: none"> <li>(→ Sample Container Requirements / Page 43): added (→ Tab. 16 Page 45) and (→ Tab. 17 Page 46).</li> </ul>

### 12.1.5 Version 09

Change No.	Changes
CHG0052462	<b>9. Other Requirements</b> <ul style="list-style-type: none"> <li>(→ Sample Container Requirements / Page 43): Added Direct Load support for Bio-Rad control tubes and Atellica TTSC, and removed VSSC.</li> </ul>

### 12.1.6 Version 08

- Added a second note to the "Delivery Information" section (→ Delivery Information / Page 8).
- Updated item #2 (0.22 micron water filters) in the "Customer Site" sub-section (→ Customer Site / Page 28) of the "Water Requirements" section.

### 12.1.7 Version 07

- Added "string and tape" and "laser level" to the list of materials (→ Materials / Page 5).
- Corrected the weights of the CH and IM analyzers and the SH in (→ Tab. 1 Page 8).
- Separated "specifications" (thermal and noise) from "requirements" (space) in the "Environmental Requirements" section (→ Environmental Requirements / Page 12).
- Added SHC thermal output to the "Thermal Output Specifications" section (→ Thermal Output Specifications / Page 12).
- Added (→ Space Requirements / Page 13).
- Changed "Is the floor acceptable?" to "Can the floor support the load of the configuration?" in step 5 of the "Environmental Requirements" section (→ Customer Site / Page 16).
- Updated the "Determining the Pitch of the Floor" section (→ Determine the Pitch of the Floor / Page 18) for clarity.

- Updated or removed the following UPS specifications:

Action	Specification	Location
Updated	Input connection	(→ Tab. 8 Page 24)
Removed	Input capacity	
Updated	Output Voltage THD	
Updated	Frequency Regulation	
Updated	Output Connections	
Updated	Recharge Time	(→ Tab. 9 Page 25)
Updated	EMC Standards	(→ Tab. 11 Page 25)

- Added a new question 5 to the “Power Requirements” section (→ Customer Site / Page 26).
- Changed deionized (DI) water to special reagent grade water (SRW) in the “Water Requirements” section (→ Water Requirements / Page 28).
- Updated text regarding 0.22 micron filter in the “Water Requirements - Customer Site” section (→ Customer Site / Page 28).
- Added a link to the “Water Quality Testing” document in the “Water Requirements” section (→ Water Requirements / Page 28).
- Added (→ Fig. 15 Page 42) to provide additional space for drawing the actual equipment location.

### 12.1.8 Version 06

- Editorial change.

### 12.1.9 Version 05

- Updated the clearance in the front of an SH in the “Space Requirements” section (→ Space Requirements / Page 13).
- Updated thermal output information (→ Thermal Output Specifications / Page 12) to include “Idle” thermal output.
- Added information regarding the LIS network settings in the “Remote Connectivity Requirements” section (→ Remote Connectivity Requirements / Page 36).

### 12.1.10 Version 04

- Added information about load per foot to the “Environmental Requirements - Customer Site” section (→ Customer Site / Page 16).

- Added the “Distances Between Feet” section (→ Distances Between Feet / Page 19).

#### 12.1.11 Version 03

- Updated the information below step 5 of the “Environmental Requirements - Customer Site” section (→ Customer Site / Page 16) to include a calculation of weight/area for each module.

#### 12.1.12 Version 02

- Initial publication to CB-DOC.

- Restricted - All documents may only be used by authorized personnel for rendering services on Siemens Healthcare Products. Any document in electronic form may be printed once. Copy and distribution of electronic documents and hardcopies is prohibited. Offenders will be liable for damages. All other rights are reserved.

[siemens-healthineers.com](http://siemens-healthineers.com)

Siemens Healthineers Headquarters  
Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen  
Germany  
Telephone: +49 9131 84-0  
[siemens-healthineers.com](http://siemens-healthineers.com)

Print No.: LDAT-000.835.01.13.02 | Replaces: LDAT-000.835.01.12.02  
Doc. Gen. Date: 02.20 | Language: English  
© Siemens Healthineers, 2020